veicomtech 802FS

veicomtech Ionic Foot Bath Detox Machine

User Manual

Model: 802FS

Introduction

The veicomtech Ionic Foot Bath Detox Machine is designed to provide a relaxing foot spa experience using advanced water electrolysis technology. This system aims to assist in detoxification and promote overall well-being by releasing ions into the water. Regular use may contribute to rejuvenation, fatigue reduction, and improved memory and sleep quality.

Important Safety Information

- Always read the instruction manual carefully before use.
- Ensure the power adapter and all connections are secure before operating the device.
- Keep the control unit away from water splashes.
- **Warm Note:** The water in the foot bath will change color due to the chemical electrolysis reaction and ions released by the product, regardless of whether feet are in the water. This is a normal operational characteristic.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult a healthcare professional for any health concerns.

Package Contents

Your veicomtech Ionic Foot Bath Detox Machine package includes the following components:

- 1 x Circular Control Unit
- 1 x Bathtub
- 1 x Power Adapter
- 2 x Ion Arrays
- 100 x Liners



Figure 1: All components included in the veicomtech Ionic Foot Bath Detox Machine package.

Setup and Operation

Follow these steps to set up and operate your ionic foot bath machine:

1. **Prepare the Tub:** Place one of the provided plastic liners inside the bathtub. This ensures easy cleanup after use.
2. **Add Water:** Pour water into the tub until it covers your ankles.
3. **Add Salt:** Add 0.4 to 0.6 ounces of salt to the water and mix thoroughly. Salt is essential for the electrolysis process.
4. **Connect Ion Array:** Connect an ion array to the circular control unit.
5. **Place Array in Water:** Carefully place the ion array into the water within the tub.
6. **Power On:** Connect the power adapter to the control unit and plug it into a power outlet. The device is designed for one-key operation, making it simple to start.
7. **Foot Placement:** Place your feet into the water for the duration of the session.
8. **Session Duration:** Foot detoxification sessions are recommended for 30 to 60 minutes per week.



Figure 2: Visual guide for the 4 steps to enjoy your foot spa.



Figure 3: The control unit features a simple one-touch start button.



Figure 4: The 30-minute timer is indicated by six lights, each representing 5 minutes of operation.

Your browser does not support the video tag.

Video 1: Official demonstration of the VEICOMTECH Foot Ionic Bath Detox Machine, showing setup and operation.

Understanding Water Color Changes

It is important to understand that the water in the foot bath will change color during operation, even without feet in the water. This is a normal phenomenon resulting from the chemical electrolysis reaction and the release of ions by the product's array. The device utilizes the principle of electrolysis, a process that breaks down water into oxygen and hydrogen. These gases are attracted by the metal array parts, leading to an oxidation reaction. The final color change is a combination of this reaction and any substances present in the water or released from the body.



Figure 5: Progression of water color change during a foot bath session, illustrating the effect of the electrolysis process.




Figure 6: The visual process of ionic detoxification, showing water changing color as the session progresses.

Maintenance and Cleaning

Maintaining your veicomtech Ionic Foot Bath is simple due to its design with disposable liners:

- **After Each Use:** Carefully remove the plastic liner from the bathtub and dispose of it. This eliminates the need to scrub the tub itself.
- **Clean Ion Array:** Rinse the ion array with clean water after each session to remove any residue. Ensure it is dry before storage.
- **Store Properly:** Store all components in a dry, safe place.

Troubleshooting

Issue	Possible Cause	Solution
Water changes color without feet in it.	Normal electrolysis reaction of the device with water and salt.	This is expected behavior and not an indication of malfunction. Continue with your session.
Device does not power on.	Power adapter not properly connected or power outlet issue.	Check all power connections. Try a different power outlet.
Ion array not bubbling or showing activity.	Insufficient salt in water or array not properly connected.	Ensure 0.4-0.6 ounces of salt are added and mixed. Verify array connection to the control unit. If using a second array, test it to rule out a faulty array.

Specifications

- **Brand:** veicomtech
- **Model Number:** 802FS
- **Product Dimensions:** 13.5 x 7 x 17.5 inches
- **Item Weight:** 1.5 Pounds
- **Material:** Plastic
- **Power Source:** Corded Electric
- **Use For:** Ankle, Foot
- **Manufacturer:** veicomtech
- **Date First Available:** August 7, 2020

Warranty and Support

veicomtech is committed to providing excellent after-sales service for your Ionic Foot Bath Detox Machine. This product is a pure water reactive foot bath system, designed for reliable use.
If you encounter any questions or issues regarding the usage of your device, please refer to this user guide or contact veicomtech customer support directly. We are dedicated to assisting you promptly and effectively.
For support, please visit the official veicomtech website or contact their customer service as indicated on the product packaging or purchase documentation.

© 2024 veicomtech. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question