BlessReach 2105

BlessReach Knee Scooter Model 2105 Instruction Manual

For Foot Injuries and Mobility Support

Overview Features Assembly Operation Maintenance & Support Troubleshooting Specifications Safety Warranty

1. Overview

The BlessReach Knee Scooter Model 2105 is designed to provide mobility support for individuals recovering from foot, ankle, or lower leg injuries, or for those who are below-the-knee amputees. It serves as an alternative to traditional crutches, offering enhanced stability and comfort during recovery.

2. Product Features

- **Adjustable Knee Rest:** Contoured and padded knee rest for comfort, designed to cradle the leg.
- **Ergonomic Handles:** Non-slip handles with an adjustable locking handbrake for safety and control.
- **Durable Wheels:** Equipped with 8-inch rubber wheels for smooth maneuvering.
- **Shock Absorption:** Features customized shock absorbers for vibration damping and noise reduction, enhancing comfort during use.
- **Foldable Design:** The rear wheel bracket and operating lever can be folded for compact storage and easy transportation.
- **Storage Bag:** Includes a detachable canvas bag with a double-layer structure for carrying personal items.
- **Weight Capacity:** Supports individuals up to 300 lbs.



Image: Customized shock absorbers for vibration damping and noise reduction, supporting up to 300 lbs.



Image: Ergonomic non-slip handle with a one-touch parking function for secure braking.



Image: 8-inch all-terrain non-slip tires suitable for various surfaces like gravel, dirt, and sand roads.



Image: Close-up of the knee pad, highlighting its breathable comfort fabric and pressure relief design.

3. Setup and Assembly

Follow these steps to assemble your BlessReach Knee Scooter Model 2105. It is recommended to review the full installation video for detailed guidance.

1. **Unfold the Frame:** Carefully unfold the main frame of the knee scooter. Ensure all locking mechanisms click into place securely.
2. **Install Front Wheels:** Attach the front wheel assembly to the main frame. Ensure all screws are tightened before moving the wrench.
3. **Adjust Handlebar Height:** Adjust the handlebar to a comfortable height, typically where your hands can rest naturally while standing upright. Secure the adjustment lever.
4. **Adjust Knee Rest Height:** Adjust the knee rest to a height that allows your injured leg to be comfortably bent at a 90-degree angle while your foot is off the ground. Secure the adjustment lever.
5. **Attach Storage Bag:** Secure the detachable canvas storage bag to the front of the scooter.

Your browser does not support the video tag.

Video: Detailed installation guide for the BlessReach Knee Scooter Model 2105.

Your browser does not support the video tag.

Video: Instructions on how to correctly unfold the knee scooter frame.

Your browser does not support the video tag.

Video: Precautions for installing the front wheel crossbar buckle.



Image: The knee scooter folded for transport, demonstrating its compact design for placement in a car.

4. Operating Instructions

To operate the BlessReach Knee Scooter, place the injured leg comfortably on the padded knee rest. Use your uninjured leg to propel the scooter forward. Steer using the handlebars. Always maintain a firm grip on the handlebars for stability.

Braking

The knee scooter is equipped with an ergonomic adjustable locking handbrake. To slow down or stop, gently squeeze the brake levers on the handlebars. To engage the parking brake, squeeze the brake lever and press the small button to lock it in place. To release, squeeze the lever again. Ensure the brakes are functioning correctly before each use.

Your browser does not support the video tag.

Video: Instructions on how to adjust the brake tension for optimal performance.

5. Maintenance

- **Regular Inspection:** Periodically check all bolts, screws, and locking mechanisms to ensure they are tight and secure.
- **Cleanliness:** Wipe down the frame and knee rest with a damp cloth as needed. Avoid harsh chemicals that may damage the finish or materials.
- **Wheel Check:** Inspect the wheels for any signs of wear or damage. Ensure they rotate freely.
- **Brake Function:** Test the brakes regularly to confirm they engage and disengage properly. Adjust brake tension if necessary (refer to operating instructions).

6. Troubleshooting

- **Squeaking Noises:** If the scooter develops squeaking noises, check all moving parts and joints. Lubrication may be required. Ensure all fasteners are tight.
- **Loose Brakes:** If brakes feel loose or ineffective, adjust the brake tension as shown in the operating instructions video.
- **Unstable Knee Pad:** If the knee pad swivels or feels insecure, check the attachment points and ensure they are properly tightened.
- **Difficulty Steering:** Ensure the handlebar is securely fastened and the front wheel assembly is correctly installed and tightened.

7. Specifications

Feature	Detail
Model Number	2105
Brand	BlessReach
Weight Capacity	300 pounds
Item Weight	22 pounds
Wheel Size	8 inches (front and rear)
Material Type	Metal
Minimum Height Recommendation	4.6 inches
Maximum Height Recommendation	6.6 inches
Package Dimensions	20.58 x 14.58 x 10 inches



Image: Diagram illustrating the key dimensions of the BlessReach Knee Scooter, including handlebar and knee rest height adjustments.

8. Safety Information

- Always read and understand the entire instruction manual before using the knee scooter.
- Do not exceed the maximum weight capacity of 300 lbs.
- Ensure all parts are securely assembled and tightened before each use.
- Use the knee scooter on stable, level surfaces. Avoid uneven terrain, stairs, or steep inclines.
- Always engage the parking brake when stationary to prevent unintended movement.
- Consult with a healthcare professional regarding appropriate use and any specific limitations related to your injury.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

9. Warranty and Support

For warranty information or product support, please refer to the documentation included with your purchase or contact BlessReach customer service directly. Keep your purchase receipt for warranty claims.

© 2026 BlessReach. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question