Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

YUYTE YUYTE

YUYTE Portable Ionic Foot Bath SPA Instruction Manual

Model: YUYTE

1. Product Overview

The YUYTE Portable Ionic Foot Bath SPA is designed to facilitate detoxification through the principle of electrolysis. It generates negative ions in water, creating an active water environment. This process aims to assist in the rearrangement of cell ions, improving cell membrane permeability for detoxification and supporting normal cellular function.
This compact and portable device is suitable for home use or travel, requiring an external foot basin or container for operation.



Figure 1.1: The YUYTE Portable Ionic Foot Bath SPA main unit and ion head, shown with feet soaking in a bowl of water with rose petals.

2. Safety Information

Important Safety Precautions:

- Not recommended for individuals with pacemakers.
- Not recommended for pregnant individuals.
- Not recommended for individuals with implanted organs.
- Not recommended for patients with epilepsy.
- Not recommended for individuals with open wounds on their feet.
- Not recommended for individuals receiving any form of radiation or chemotherapy.
- Not recommended for patients with type 1 diabetes.

Always consult with a healthcare professional before using this device, especially if you have pre-existing medical conditions.

3. Package Contents

Please verify that all items listed below are included in your package:

- 1 x Main Machine (Ionic Detox Machine)
- 1 x Ion Head (Array)
- 1 x Power Adapter
- 1 x Spoon (for salt measurement)
- 1 x User Manual (this document)
- 1 x Charger (likely referring to the power adapter)



Figure 3.1: All components included in the YUYTE Portable Ionic Foot Bath SPA package.

4. Setup Instructions

1. **Prepare the Basin:** Obtain a suitable foot basin or any container large enough for soaking your feet. The product does not include a basin.
2. **Place the Main Machine:** Position the main machine on a stable, dry surface near the basin.
3. **Connect the Ion Head:** Connect the ion head (array) cable to the "ARRAY" port on the main machine.
4. **Place Ion Head in Water:** Submerge the ion head completely in the water within your chosen foot basin.
5. **Add Salt:** Add 1 to 3 grams of salt to the water. The salt helps to increase the conductivity of the water, which is

essential for the electrolysis process.

6. **Connect Power:** Connect the power adapter to the "DC INPUT" port on the main machine. Plug the adapter into a standard electrical outlet. Ensure the voltage is compatible with your region (110-240V).
7. **Power On:** Press the "On / Off" button on the main machine to power on the device.



Figure 4.1: The main unit and ion head connected, illustrating the setup.

5. Operating Instructions

1. **Initiate Session:** After completing the setup, ensure your feet are clean and comfortably placed in the basin with the ion head.
2. **Set Time:** The device allows you to set the session time according to your needs. Refer to the device's interface for time adjustment options.
3. **Monitor Salt Concentration:** The device's design allows for easy observation of salt concentration changes in the water, which may indicate the progress of the session.
4. **Session Duration:** For optimal results, a typical session duration is recommended. The manufacturer suggests using the device twice a week. After approximately 12 sessions, users may begin to notice significant effects.
5. **End Session:** Once the desired session time is complete, or if you wish to stop early, press the "On / Off" button to power off the device.
6. **Remove Feet:** Carefully remove your feet from the basin.

The ionic detox process aims to promote the body's natural detoxification mechanisms and enhance overall well-being by

supporting cellular function and energy levels.



Figure 5.1: The device is designed to promote metabolism, improve blood circulation, enhance healthy cell regeneration, and strengthen the immune system.



Figure 5.2: Potential benefits include increased energy, reduced stress, pain relief, improved memory and sleep, and reduction of certain skin conditions.

6. Maintenance

Proper maintenance ensures the longevity and effective operation of your Ionic Foot Bath SPA.

- **Cleaning the Ion Head:** After each use, rinse the ion head thoroughly with clean water to remove any residue. Ensure it is completely dry before storage.
- **Cleaning the Main Machine:** The main machine has a hard plastic casing. Wipe it with a soft, damp cloth. Do not submerge the main machine in water or use abrasive cleaners.
- **Storage:** Store the device and its components in a cool, dry place away from direct sunlight and extreme temperatures.

7. Troubleshooting

If you encounter issues with your YUYTE Portable Ionic Foot Bath SPA, refer to the following common troubleshooting

steps:

- **Device Not Powering On:**
  - Ensure the power adapter is securely connected to both the device and a working electrical outlet.
  - Verify that the electrical outlet is functional by plugging in another appliance.
  - Check the power adapter for any visible damage.
- **No Bubbling/Activity in Water:**
  - Confirm that the ion head is fully submerged in the water.
  - Ensure that 1 to 3 grams of salt have been added to the water. Insufficient salt will prevent proper conductivity.
  - Check that the ion head cable is securely connected to the "ARRAY" port.
- **Unusual Odor or Noise:**
  - If you notice any unusual odors or loud noises, immediately power off the device and unplug it. Contact customer support if the issue persists.

For issues not covered here, please contact YUYTE customer support.

8. Specifications

Feature	Specification
Item Type	Ionic Detox Machine
Material	ABS Plastic
Voltage	110-240V
Main Machine Size	Approx. 12 x 12 x 5.5 cm (4.7 x 4.7 x 2.2 inches)
Weight	1.57 Pounds (approx. 0.71 kg)
Color	White
Power Source	Corded Electric
Model Number	YUYTE
Manufacturer	Cimenn



Figure 8.1: Dimensions of the main unit.

9. Warranty and Support

For product support, inquiries, or warranty information, please visit the official YUYTE store or contact their customer service directly.

YUYTE Store: https://www.amazon.com/stores/YUYTE/page/2AE18F32-DFA7-41C8-8423-0F33D3A18AE4

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

© 2024 YUYTE. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question

and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0987NMBXQ

More documents and user manuals: **https://manuals.plus/**