Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

waxkiss PWH-5000B

Waxkiss PWH-5000B Paraffin Wax Machine User Manual

Model: PWH-5000B

1. Important Safety Instructions

Please read all instructions carefully before using the Waxkiss PWH-5000B Paraffin Wax Machine. Keep this manual for future reference.

- **Avoid contact with eyes.** If wax comes into contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Do not immerse the machine in water or any other liquid.
- Always unplug the machine when not in use, before cleaning, or if a malfunction occurs.
- Do not operate the machine with a damaged cord or plug.
- Keep out of reach of children.
- Ensure the machine is placed on a stable, heat-resistant surface.
- Do not use if you have open wounds, skin infections, or severe skin conditions. Consult a physician if you have concerns.
- Always test the wax temperature on a small area of your skin before full immersion to prevent burns.

2. Package Contents

Verify that all items are present in the package:

- 1 x Waxkiss Paraffin Wax Machine (5000ml capacity)
- 4 x Packages of Lavender Paraffin Wax
- 1 x Pair of Cotton Mitts for Hands
- 1 x Pair of Cotton Booties for Feet
- Disposable PE Mitts for Hands (quantity may vary)
- Disposable PE Booties for Feet (quantity may vary)
- 1 x Heat Insulating Grid
- 1 x User Manual (this document)



Image 2.1: Contents of the Waxkiss Paraffin Wax Machine kit, including the machine, paraffin wax blocks, cotton mitts, cotton booties, and heat insulating grid.

3. Product Overview

The Waxkiss PWH-5000B is a large-capacity paraffin bath designed for professional spa treatments and arthritis relief at home. It features precise digital temperature control and a spacious design to accommodate hands, feet, and elbows.

Key Features:

- **Large Capacity:** 5000ml pot, suitable for hands, feet, and elbows.
- **Precision Temperature Control:** Digital display with a temperature range of 86°F to 185°F (30°C to 85°C).
- **Timer Function:** Adjustable time setting from 1 to 99 minutes.
- **User-Friendly Interface:** Digital display and control buttons for easy operation.
- **Durable Construction:** Made from safe materials for reliable use.



Image 3.1: Top-down and side view of the Waxkiss Paraffin Wax Machine illustrating its dimensions, including a 13-inch length and 7-inch width for the opening, and an overall height of 8 inches.

4. Setup

1. **Unpack:** Carefully remove all components from the packaging.
2. **Placement:** Place the paraffin wax machine on a flat, stable, and heat-resistant surface. Ensure there is adequate ventilation around the unit.
3. **Clean:** Before first use, wipe the inside of the wax pot with a damp cloth to remove any dust. Ensure it is completely dry before proceeding.
4. **Add Wax:** Place the paraffin wax blocks directly into the machine's pot. Do not add water or any other liquids.
5. **Connect Power:** Plug the power cord into a standard electrical outlet.

5. Operating Instructions

5.1 Melting the Wax

1. **Power On:** Press the 'ON/OFF' button to turn on the machine. The digital display will illuminate.
2. **Select Mode:** The machine has two working modes: Temperature Heating Mode and Timed Heating Mode. Press the 'SET' button to toggle between these modes.
3. **Set Temperature (Temperature Heating Mode):** If in Temperature Heating Mode, use the '+' or '-' buttons to adjust the desired temperature. For melting wax, a higher temperature like 185°F (85°C) is recommended initially.
4. **Set Time (Timed Heating Mode):** If in Timed Heating Mode, use the '+' or '-' buttons to set the heating duration (1-99 minutes). The machine will automatically maintain constant heating according to the wax amount.
5. **Melt Wax:** Allow the wax to melt completely. This typically takes about 50 minutes for 2 lbs of wax at 185°F (85°C). The machine will emit three 'Di, Di, Di' sounds when the wax reaches 122°F (50°C), indicating it is ready for use.
6. **Insert Grid:** Once the wax is fully melted, place the heat insulating grid at the bottom of the pot. This prevents direct contact with the heating element.



Image 5.1: Diagram illustrating the two working modes of the Waxkiss Paraffin Wax Machine: Temperature Heating Mode and Timed Heating Mode, with controls for setting temperature and time.

5.2 Paraffin Treatment Application

1. **Test Temperature:** Before immersing your hand or foot, carefully test the wax temperature on a small, less sensitive area of your skin (e.g., wrist) to ensure it is comfortable and not too hot. The ideal temperature for application is typically around 125-130°F (52-54°C). Adjust using the '+' or '-' buttons if needed.
2. **Prepare Skin:** Ensure your hands, feet, or elbows are clean and dry.
3. **Dip:** Slowly dip your hand, foot, or elbow into the melted paraffin wax. Remove it for a few seconds to allow the wax to cool slightly and solidify, then dip again. Repeat this process 3-5 times to build up a thick layer of wax.
4. **Cover:** Immediately after the final dip, cover the waxed area with a disposable PE mitt/bootie.
5. **Insulate:** Place the cotton mitt or booties over the disposable cover to retain heat.
6. **Relax:** Allow the wax to remain on your skin for 10-15 minutes. This allows the heat and moisturizing properties to penetrate.
7. **Remove Wax:** Gently peel off the solidified paraffin wax. It should come off easily, taking dead skin cells with it. Discard the used wax.
8. **Moisturize:** After removing the wax, you may apply a moisturizer to further hydrate your skin.




Image 5.2: A six-step visual guide demonstrating the paraffin wax treatment process: 1. Put wax into pot, 2. Dip hand after melting and placing grid, 3. Wrap hand with glove liner, 4. Put on cotton mitt, 5. Remove wax after 15 minutes, 6. Resulting soft skin.

5.3 Power Off

When finished, press the 'ON/OFF' button to turn off the machine and unplug it from the power outlet. Allow the machine and remaining wax to cool completely before handling or storing.

6. Maintenance and Cleaning

1. **Unplug and Cool:** Always ensure the machine is unplugged and completely cool before cleaning.
2. **Remove Wax:** If you need to clean the pot, allow the wax to solidify. Once solid, you can often gently lift the entire block of wax out. If wax remains, you can reheat it slightly to soften and then wipe it out with a paper towel or cloth.
3. **Clean Interior:** Use a soft, damp cloth to wipe the interior of the wax pot. Do not use abrasive cleaners or solvents. Ensure the pot is completely dry before next use.
4. **Clean Exterior:** Wipe the exterior of the machine with a soft, damp cloth.

5. **Storage:** Store the machine in a cool, dry place, away from direct sunlight and moisture.

7. Troubleshooting

- **Machine Does Not Turn On:**
  - Check if the power cord is securely plugged into both the machine and the electrical outlet.
  - Ensure the power outlet is functional by testing it with another appliance.
- **Wax Not Melting:**
  - Verify that the machine is powered on and the temperature is set appropriately (e.g., 185°F/85°C for initial melting).
  - Allow sufficient time for the wax to melt, especially if using a large quantity.
- **Wax Too Hot/Cold:**
  - Adjust the temperature setting using the '+' or '-' buttons. Allow a few minutes for the wax temperature to stabilize after adjustment.
  - Always test the wax temperature before application.

If you encounter issues not listed here, please contact customer support for assistance.

8. Specifications

Model Number	PWH-5000B
Brand	waxkiss
Capacity	5000ml
Temperature Range	86°F - 185°F (30°C - 85°C)
Timer Range	1 - 99 minutes
Package Dimensions	16.81 x 11.54 x 7.91 inches
Item Weight	6.75 pounds
UPC/GTIN	734257215466

9. Warranty and Support

For warranty information and customer support, please refer to the documentation included with your purchase or visit the official Waxkiss website. If you have any questions or require assistance with your Waxkiss PWH-5000B Paraffin Wax Machine, please contact the manufacturer's customer service department.

You can typically find support contact details on the product packaging or by visiting the Waxkiss Store on Amazon.

© 2023 Waxkiss. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B099RD8LSQ

More documents and user manuals: **https://manuals.plus/**