LeCuag 805 D

Ionic Foot Bath Detox Machine User Manual

Brand: LeCuag | Model: 805 D

Introduction

This user manual provides comprehensive instructions for the safe and effective operation and maintenance of your LeCuag Ionic Foot Bath Detox Machine. Please read this manual thoroughly before first use to ensure proper functionality and to maximize the benefits of your device.

Important Safety Information

WARNING: Do not use this appliance if your feet are swollen or inflamed or if you have any skin irritations. Always consult with a healthcare professional before using this device, especially if you have pre-existing medical conditions, are pregnant, or have electronic implants.
Ensure the machine's voltage is 110V before use. Always check that the machine is properly packaged upon receipt.

Package Contents

Please verify that all items listed below are included in your package:

- 1 x LeCuag Ionic Foot Bath Detox Machine (Main Unit)
- 2 x Ion Arrays
- 2 x Wrist Straps
- 1 x Far Infrared Belt
- 1 x Power Cord
- 2 x Electrode Patches
- User Manual (this document)
- Aluminum Alloy Carrying Case



Image: The complete LeCuag Ionic Foot Bath Detox Machine package, including the main unit, ion arrays, wrist straps, infrared belt, and other accessories, neatly arranged in its aluminum alloy carrying case.

Product Features

- **Personalized Detoxification:** Features 8 adjustable foot detox modes to cater to different physiological needs and desired intensity.
- **Flexible Session Times:** Offers humanized time settings from 10 to 30 minutes, allowing for efficient use.
- **Integrated Massage Modes:** Includes 5 massage modes for enhanced relaxation and comfort during sessions.
- **Smart Interruption Function:** Allows users to pause and resume detox sessions at any time.
- **Safety Auto-Stop:** The machine will beep and stop operation if the amperage exceeds 2.5A, indicating a need to adjust water or replace the array.
- **Dual User Capability:** Designed to accommodate 1-2 persons simultaneously, promoting shared wellness experiences.
- **Far Infrared Belt:** Helps improve circulation, relieve coldness in hands and feet, and may assist with menstrual pain, constipation, and metabolism.
- **Clear Digital Display:** Equipped with a large colored screen that clearly shows mode, time, and current amperage, making it

user-friendly, especially for the elderly.

- **Portable Design:** Comes with a luxury Aluminum alloy box for easy transport and travel.



Image: Visual representation of the 3-in-1 benefits of the LeCuag Ionic Foot Bath Detox Machine, highlighting its ability to reduce stress, relieve pain, improve sleep, remove heavy metals, and aid in weight management.



Image: A user relaxing while utilizing the Far Infrared Belt, illustrating its benefits for circulation and overall well-being.



Image: A clear diagram explaining the 8 distinct detoxification modes, detailing their positive and negative ion percentages and recommended usage for different body types or conditions.

Setup Instructions

Follow these steps to set up your Ionic Foot Bath Detox Machine for use:

1. **Prepare the Basin:** Place a suitable foot bath basin (not included) on a flat, stable surface. Line the basin with a plastic liner for easier cleaning, if desired.
2. **Add Water:** Fill the basin with warm water, ensuring it covers your feet comfortably.
3. **Add Salt (Optional but Recommended):** Add a small amount of Himalayan salt to the water. This helps improve conductivity and the detoxification process. Stir until dissolved.
4. **Connect the Ion Array:** Place one or two ion arrays into the water in the basin. Connect the array cable(s) to the designated "ARRAY" port(s) on the main unit.
5. **Wear the Wristband:** Secure the wristband(s) firmly around your wrist(s). Ensure the metal sheet of the wristband makes good contact with your skin. Connect the wristband cable(s) to the designated "WRIST BAND" port(s) on the main unit.

6. **Connect Far Infrared Belt (Optional):** If using, wrap the Far Infrared Belt around the desired area (e.g., waist, back) and connect its cable to the "LF" (Low Frequency) port on the main unit.
7. **Connect Power:** Plug the power cord into the main unit and then into a standard 110V electrical outlet.
8. **Power On:** Press the power button on the main unit to turn it on. The display should illuminate. If the display is dark after 10 seconds, press "Start/Stop" to wake it up.



Image: Two users enjoying a foot bath session, demonstrating the dual-user capability of the machine. Note that the foot basins are not included.

Operating Instructions

Once the machine is set up, follow these steps to begin your detox session:

1. **Place Your Feet:** Gently place your feet into the water in the basin. Ensure your feet are fully submerged and the wristband is securely worn. The machine will not operate if feet are not in the water or the wristband is not properly connected.
2. **Select Detox Mode:** Use the "MODE" button to cycle through the 8 available detox modes. Refer to the "8 Modes As Following" diagram for details on each mode's ion concentration.

3. **Set Time:** Use the "TIME" buttons (+/-) to adjust the session duration (10-30 minutes).
4. **Adjust Far Infrared Belt (Optional):** If using the Far Infrared Belt, use the "A/B Belt" buttons to select the desired intensity (L, M, H).
5. **Start Session:** Press the "START/STOP" button to begin the detoxification process. The display will show the active mode, time remaining, and current amperage.
6. **Monitor Display:** Observe the display for voltage and current readings. If the amperage exceeds 2.5A, the machine will beep and stop. This indicates the salt concentration might be too high or the array needs cleaning/replacement.
7. **Massage Function (Optional):** Press the "ON/OFF" button for the massage function. Use the "+/-" buttons to adjust the massage intensity.
8. **End Session:** The machine will automatically stop when the set time expires. You can also press "START/STOP" to pause or end the session manually.
9. **After Session:** Remove your feet from the basin, unplug the machine, and proceed with cleaning as described in the Maintenance section.



Image: An annotated diagram of the control panel, explaining the function of each button (START/STOP, GRADE A/B MODE, A/B TIME, A/B Belt, ON/OFF, +/-) and the information displayed on the screen, including voltage and ion concentration.

Maintenance

Proper maintenance ensures the longevity and optimal performance of your LeCuag Ionic Foot Bath Detox Machine:

- **After Each Use:**
  - Unplug the machine from the power outlet.
  - Remove the ion array(s) from the basin.
  - Discard the water and any plastic liners from the basin.
  - Rinse the ion array(s) thoroughly under running water to remove any residue. You may use a soft brush to gently clean the coils.
  - Wipe down the main unit with a soft, damp cloth. Do not immerse the main unit in water.
  - Ensure all components are completely dry before storing.
- **Storage:** Store the machine and its accessories in the provided aluminum alloy case in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Array Replacement:** Ion arrays are consumable parts and will degrade over time. Replace them when you notice a significant decrease in bubbling or if the machine frequently stops due to high amperage even with proper water and salt levels.

Troubleshooting

If you encounter any issues with your LeCuag Ionic Foot Bath Detox Machine, please refer to the common problems and solutions below:

Problem	Possible Cause	Solution
Machine does not turn on or display is dark.	No power, machine in sleep mode, or loose connection.	Ensure power cord is securely plugged into the machine and outlet. Press "Start/Stop" to wake up the display. Verify voltage is 110V.
Voltage display shows "0" and/or wristband/foot bath icons flash.	Wristband not worn or loose connection; feet not in water.	Ensure the wristband is worn securely with the metal sheet touching your wrist. Place both feet fully into the water.
Machine beeps and stops (Amperage exceeds 2.5A or 3.0A).	Salt concentration too high, dirty array, or worn-out array.	Change the water and use pure water. Clean the ion array thoroughly. If the issue persists, replace the ion array.
Ion array does not bubble sufficiently.	Insufficient salt, dirty array, or worn-out array.	Add a small amount of salt to the water. Clean the ion array. If still no bubbles, replace the ion array.
Tingling sensation on wrist when wearing wristband.	Dry skin contact with metal sheet.	Stop the machine. Wet your wrist or apply a small amount of hand cream to the area where the metal sheet touches. You can also place a piece of wet paper under the metal sheet.

Your browser does not support the video tag.

Video: A troubleshooting guide for the LeCuag Dual Ionic Detox Machine, covering common issues such as the machine not starting, voltage display errors, beeping due to high amperage, and lack of array bubbling. It also demonstrates how to properly wear the wristband and clean the array.

Specifications

Attribute	Value

Attribute	Value
Product Dimensions	17 x 13 x 7 inches
Item Weight	16.07 Pounds
Manufacturer	LeCuag
Model Number	805 D (Inferred from video)
ASIN	B09R6RZR2S
UPC/GTIN	735209275576

Warranty and Support

For warranty information, technical support, or any inquiries regarding your LeCuag Ionic Foot Bath Detox Machine, please contact the manufacturer directly through their official website or the retailer where the product was purchased. Keep your purchase receipt as proof of purchase for warranty claims.

© 2024 LeCuag. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question