Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

onerbl 6-Volt AC Adapter

onerbl 6-Volt AC Adapter Instruction Manual

For CVS Automatic Blood Pressure Monitors (Models 306409, 271245, 271244, 344532)

1. Introduction

This manual provides essential information for the safe and effective use of your onerbl 6-Volt AC Adapter. This adapter is designed as a replacement power supply for various CVS Automatic Blood Pressure Monitor models, including 306409, 271245, 271244, and 344532. It features a universal input voltage range (AC 100-240V) and incorporates safety protections against over voltage and overheating.



Image 1.1: Front and side view of the onerbl 6-Volt AC Adapter, showing the two-prong plug and the barrel connector cable.

2. Setup

1. **Identify Device Power Input:** Locate the power input port on your CVS Automatic Blood Pressure Monitor. Ensure it requires a 6-volt AC adapter with a barrel connector.
2. **Connect Adapter to Device:** Insert the barrel connector end of the onerbl AC adapter firmly into the power input port of your blood pressure monitor.
3. **Connect to Power Outlet:** Plug the two-prong end of the AC adapter into a standard AC 100-240V wall outlet. Ensure the outlet is easily accessible and in good condition.




Image 2.1: The onerbl AC Adapter plugged into a power strip, demonstrating a typical connection scenario.



Image 2.2: The onerbl AC Adapter plugged directly into a wall outlet.

3. Operating Instructions

Once the onerbl 6-Volt AC Adapter is correctly connected to both your blood pressure monitor and a power outlet, your device should receive power. The adapter operates automatically, converting the AC wall power to the required 6-volt DC power for your monitor. No further user interaction is required for the adapter itself during operation.

4. Maintenance

- **Cleaning:** Keep the adapter clean and free from dust. Use a dry, soft cloth to wipe the surface. Do not use liquid cleaners.
- **Storage:** Store the adapter in a cool, dry place when not in use. Avoid extreme temperatures and high humidity.
- **Cable Inspection:** Regularly inspect the power cord and connector for any signs of damage, such as fraying, cuts, or bent pins. Discontinue use if any damage is observed.
- **Ventilation:** Ensure the adapter has adequate ventilation during use. Do not cover it or place it in enclosed spaces that could lead to overheating.

5. Troubleshooting

- **Device Not Powering On:**
  - Ensure the adapter's barrel connector is fully inserted into the blood pressure monitor.
  - Verify the adapter is securely plugged into a working wall outlet. Try a different outlet if necessary.
  - Confirm that your blood pressure monitor model is listed as compatible with this 6-volt adapter.
- **Adapter Overheating:**
  - Disconnect the adapter from the power source and the device immediately.
  - Check for any obstructions around the adapter that might impede airflow.
  - If overheating persists after ensuring proper ventilation, discontinue use and contact support.
- **Intermittent Power:**
  - Check all connections for looseness.
  - Inspect the cable for any visible damage or kinks.

The adapter is equipped with built-in protections against over voltage and over heat to ensure safe operation.



Image 5.1: The onerbl AC Adapter features multiple safety protections, including short-circuit, over-voltage, and over-current protection.

6. Specifications

Feature	Detail
Brand	onerbl
Input Voltage	AC 100-240V, 50/60Hz
Output Voltage	DC 6V (specific current not provided, but implied by compatibility)

Connector Type	Barrel Connector
Compatible Devices	CVS Automatic Blood Pressure Monitors (Models 306409, 271245, 271244, 344532)
Cable Length	Approximately 76 inches (6.33 feet)
Item Weight	8 Ounces
Main Power Connector Type	2 Pin (US Standard)
Special Features	Travel-friendly, Over Voltage Protection, Over Heat Protection
Manufacturer	Shenzhen Waweis Technology Co., Ltd



Image 6.1: Dimensions of the onerbl AC Adapter and its 76-inch cable length.



Image 6.2: Illustration of the adapter's durable design, featuring copper wire and a PVC jacket.

7. Support

For further assistance, technical support, or warranty inquiries regarding your onerbl 6-Volt AC Adapter, please contact the manufacturer or your retailer. Refer to your purchase documentation for specific contact details.

© 2024 onerbl. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0B152FVZ6

More documents and user manuals: **https://manuals.plus/**