Paingone B0B4XPRJTF

Paingone Easy Replacement Gel Pads Instruction Manual

Model: B0B4XPRJTF

1. Product Overview

The Paingone Easy Replacement Gel Pads are designed exclusively for use with the Paingone Easy device (sold separately). These premium-quality, long-life adhesive electrode pads are engineered to deliver stimulating TENS (Transcutaneous Electrical Nerve Stimulation) therapy to various parts of the body, aiding in pain relief.
Each pack contains two replacement pads, offering up to 3 months of use based on typical application (1-2 applications per day, 15 minutes per application). Their slim and flexible design allows for discreet wear under clothing and conforms to body contours, ensuring effective therapy delivery.


Easy
This pack contains:
2x replacement electrode pads
for the Paingone Easy device
2x replacement pads
for up to 3 months use
paingone
Easy
paingone

Image 1.1: Paingone Easy Replacement Gel Pads packaging and a single pad. The packaging indicates '2x replacement pads for up to 3 months use'.



Image 1.2: A woman demonstrating the wearable nature of the Paingone Easy TENS pad on her upper back. The image highlights that the two replacement TENS pads can last up to three months of daily use with the Paingone Easy device.

2. Setup and Application

Before applying the replacement gel pads, ensure your skin is clean and dry. These pads are designed to be used exclusively with the Paingone Easy device.

1. **Attach to Device:** Connect the replacement gel pads securely to your Paingone Easy device. Ensure the connection points are firm.
2. **Apply to Skin:** Once attached to the device, firmly press the entire surface of the pads against the desired area of your skin. The flexible wings of the pads are designed to contour to the body, helping to hold them firmly in place.
3. **Placement:** Place the pads directly on the area where you experience pain. The pads are suitable for larger areas of the body.



Image 2.1: Visual guide on how to use the Paingone Easy pads. It illustrates applying the pads to the body, dampening them if stickiness is lost, and replacing them when they no longer adhere.



Image 2.2: A man placing the Paingone Easy device with pads directly on his lower back, indicating direct application to the painful area. Note: Easy TENS device not included, replacement pads only.

3. Operating Instructions

Once the Paingone Easy Replacement Gel Pads are securely attached to your Paingone Easy device and applied to your skin, operate the device according to its specific instructions. The pads facilitate the delivery of the TENS therapy.

- **Therapy Duration:** A typical application lasts 15 minutes.
- **Discreet Use:** The slim design allows the pads to be worn discreetly under clothing, enabling use at home, work, during activity, or while relaxing.
- **Targeted Relief:** The pads are designed to deliver stimulating TENS therapy directly where you need it most.




Image 3.1: A silhouette of a human body highlighting various areas where Paingone Easy TENS Pads can be placed for quick relief from aches and pains, including arms, shoulders, and legs. It mentions adjustable intensity and modes, and 15-minute therapy.



Image 3.2: A man holding his neck, illustrating the customizable relief offered by the Paingone Easy system. It highlights 12 adjustable intensity levels, a 15-minute automatic therapy program, and 3 different therapy modes (soft pulses, strong pulses, or both).

4. Maintenance

Proper maintenance of your Paingone Easy Replacement Gel Pads will help extend their lifespan and ensure effective therapy.

- **Restoring Stickiness:** If the pads begin to lose their adhesive quality, they can often be refreshed. Lightly dampen the adhesive surface of the pads with a small amount of water. Allow them to air dry slightly until they regain tackiness.
- **Replacement:** When the pads no longer adhere to the skin even after dampening, or if they show signs of wear and tear, they should be replaced with new Paingone Easy Replacement Gel Pads.
- **Storage:** Store the pads in a cool, dry place, ideally re-applying the protective film or placing them back in their original

packaging when not in use to prevent drying out and maintain hygiene.

5. Troubleshooting

If you encounter issues with your Paingone Easy Replacement Gel Pads, consider the following:

- **Pads Not Sticking:** Ensure your skin is clean and free of lotions, oils, or sweat. Try dampening the pads as described in the Maintenance section. If they still do not adhere, they may need to be replaced.
- **No Sensation/Weak Sensation:** Verify that the pads are securely attached to the Paingone Easy device and that the entire surface of the pads is in full contact with your skin. Ensure the Paingone Easy device itself is functioning correctly and charged.
- **Skin Irritation:** Discontinue use immediately if you experience any skin irritation, redness, or discomfort. Ensure the pads are clean before each use. Do not apply to broken or irritated skin.

For further assistance, please refer to the Paingone Easy device manual or contact Paingone customer support.

6. Specifications

Attribute	Detail
Product Name	Paingone Easy Replacement Gel Pads
Model Number	B0B4XPRJTF
Quantity	2 Count
Estimated Lifespan	Up to 3 months (based on 1-2 applications/day, 15 mins/application)
Product Dimensions	5.9 x 2.9 x 0.75 inches
Item Weight	1.13 ounces
Manufacturer	Paingone
Date First Available	September 9, 2022

7. Important Information

Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

8. Support and Contact

For product support, warranty information, or to purchase additional replacement pads, please visit the official Paingone website or contact their customer service department.
You can find more information and contact details by visiting thePaingone Store on Amazon.

© 2024 Paingone. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question