LeCuag L-056

LeCuag Ionic Foot Bath Spa Machine L-056 User Manual

1. Introduction

This manual provides essential instructions for the safe and effective use of your LeCuag Ionic Foot Bath Spa Machine, Model L-056. Please read this manual thoroughly before operating the device and retain it for future reference.

2. Safety Information

Adherence to the following safety guidelines is crucial for preventing injury and damage to the device:

- Use the device strictly according to the manufacturer's instructions.
- This device is not intended for use by children or individuals with certain medical conditions. Consult a healthcare professional if you have concerns.
- Ensure all electrical connections are secure and the power source meets the device's requirements to prevent electrical hazards.
- Do not immerse the control unit in water.
- Unplug the device from the power outlet when not in use and before cleaning.
- Keep the device away from heat sources and direct sunlight.
- Do not operate the device if it has a damaged cord or plug, or if it is not working properly.

3. Package Contents

Verify that all items are present in the package:

- 1 x Foot Basin
- 1 x Controller Unit
- 2 x Ion Arrays
- 100 x Foot Bath Liners
- 1 x Power Adapter
- 1 x User Manual (this document)



Image: The LeCuag Ionic Foot Bath Spa Machine, including the foot basin, controller unit, two ion arrays, power adapter, and a pack of foot bath liners.

4. Product Overview

The LeCuag Ionic Foot Bath Spa Machine is designed for personal foot care. It consists of a foot basin, a control unit, and ion arrays that facilitate an electrochemical reaction in water.

Key Components:

- **Foot Basin:** Holds water and the user's feet during operation.
- **Controller Unit:** Manages the electrical current and session timing. Features a simple on/off button and indicator lights.
- **Ion Arrays:** Submerged in water, these components facilitate the ionic process.
- **Foot Bath Liners:** Disposable liners for hygienic use.

5. Setup

Follow these steps to set up your LeCuag Ionic Foot Bath Spa Machine:

1. Place the foot basin on a flat, stable surface.
2. Insert a disposable foot bath liner into the basin.
3. Fill the basin with warm tap water, ensuring the water level covers the ion array when submerged. Do not use purified or distilled water, as minerals are necessary for the device's function.
4. Connect the ion array(s) to the controller unit.

5. Place the ion array(s) into the water inside the basin.
6. Connect the power adapter to the controller unit, then plug the adapter into a standard electrical outlet.



Image: A visual guide demonstrating the setup process, including connecting the power, placing the ion array, adding water, and inserting feet into the basin.

6. Operation

The device features a simplified operation with a pre-programmed session time.

1. After setup, ensure your feet are clean and comfortably placed in the water.
2. Press the On/Off button on the controller unit to start the session. The device has a pre-programmed 30-minute setting.
3. The device will automatically turn off after 30 minutes.
4. For optimal results, it is recommended to use the device two to three times per week.



Image: The controller unit highlighting its one-button start and 30-minute preset timer for ease of use.

Understanding Water Color Change:

The device operates by introducing an electrical current into the water via an ion array. This process causes an electrochemical reaction with the water and any dissolved minerals or salts present. This interaction results in a change in water color, which can occur even without feet submerged, due to the interaction with the metal components of the array. The water color may further darken when feet are immersed.



Image: A visual representation showing the water changing color over 5, 15, and 30 minutes, illustrating the electrochemical reaction.

7. Maintenance and Care

Proper maintenance ensures the longevity and performance of your device.

Cleaning the Ion Array:

After each use, it is important to clean the ion array. Simply spray it with water and gently brush away any residue. Regular cleaning can extend the lifespan of the ion array to approximately 40 to 50 uses.



Image: A close-up view of an ion array being cleaned with water and a brush, demonstrating the recommended maintenance procedure.

Cleaning the Foot Basin:

After each session, dispose of the used foot bath liner and rinse the foot basin with clean water. Wipe it dry with a soft cloth.

Storage:

Store the device in a cool, dry place away from direct sunlight and out of reach of children.

8. Troubleshooting

If you encounter issues with your device, refer to the following common problems and solutions:

- **Device does not power on:**
  - Ensure the power adapter is securely connected to both the controller unit and a working electrical outlet.
  - Check if the power outlet is functional by plugging in another device.

- **Water does not change color:**
  - Verify that you are using tap water, not purified or distilled water, as minerals are essential for the electrochemical reaction.
  - Ensure the ion array is properly connected and fully submerged in water.
  - The ion array may be worn out. Replace it if it has been used approximately 40-50 times.
- **Unusual noise or odor:**
  - Immediately discontinue use and unplug the device. Contact customer support for assistance.

9. Specifications

Feature	Specification
Model Number	L-056
Product Dimensions	17.48 x 13.5 x 6.57 inches
Item Weight	6.46 Pounds
Power Source	Corded Electric
Material	Metal (Ion Array), Plastic (Basin)
Color	White
Manufacturer	LeCuag

10. Warranty and Support

For warranty information or technical support, please refer to the contact details provided with your purchase documentation or visit the official LeCuag website. Keep your purchase receipt as proof of purchase.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question