Pawfan B0BWKCV2WL

Pawfan LED Facial Mask User Manual

Model: B0BWKCV2WL

Introduction

This manual provides comprehensive instructions for the safe and effective use of your Pawfan LED Facial Mask. Please read this manual thoroughly before first use and retain it for future reference.



Image: The Pawfan LED Facial Mask, a black mask with circular patterns, shown alongside its remote control. A woman is depicted wearing the mask with a red glow, indicating its light therapy function.

Important Safety Information

- Consult a healthcare professional before use if you have any pre-existing skin conditions, light sensitivity, or are taking medications that cause photosensitivity.
- Do not use on broken or irritated skin.
- Keep out of reach of children.
- Avoid direct eye exposure to the LED lights.
- Do not immerse the device in water.
- Use only the provided power adapter and remote control.
- Discontinue use if you experience any discomfort or adverse reactions.

Package Contents

Upon opening your Pawfan LED Facial Mask package, please ensure all items are present:

- LED Facial Mask
- Remote Control
- Power Adapter
- User Manual (this document)

Setup Instructions

1. **Unpack:** Carefully remove all components from the packaging.
2. **Connect Power:** Connect the power adapter to the remote control, and then connect the remote control to the LED Facial Mask.
3. **Plug In:** Plug the power adapter into a suitable electrical outlet.
4. **Adjust Straps:** Place the mask on your face and adjust the head straps for a comfortable and secure fit. Ensure the mask is not too tight.



Image: The Pawfan LED Facial Mask being held by a person, with the remote control visible. The mask is illuminated with both red and blue light, indicating different treatment modes.

Operating Instructions

The Pawfan LED Facial Mask is an over-the-counter device. The red light mode is intended for the treatment of periorbital wrinkles.

1. **Power On:** Press the power button on the remote control to turn on the mask.
2. **Select Light Mode:** Use the mode selection button on the remote to cycle through available light colors (e.g., Red Light, Blue Light).
3. **Set Treatment Time:** Use the timer buttons (e.g., 5min, 10min, 15min) on the remote control to select your desired treatment duration.
4. **Begin Treatment:** Relax while wearing the mask for the selected duration.
5. **Power Off:** The mask will automatically turn off after the set time, or you can manually turn it off by pressing the power button.

Recommended Usage:

- For optimal results, use consistently for 10-15 minutes daily.
- Ensure your face is clean and dry before each session.



Image: A collage showing the Pawfan LED Facial Mask in use during various activities, including reading, sleeping, and exercising, highlighting its lightweight and portable design.

Care and Maintenance

- **Cleaning:** Gently wipe the mask with a soft, damp cloth after each use. Do not use abrasive cleaners or harsh chemicals. Ensure the mask is completely dry before storing.
- **Storage:** Store the mask in a cool, dry place away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.
- **Handling:** Handle the mask and remote control with care to avoid dropping or physical impact.

Troubleshooting

Problem	Possible Cause	Solution
Mask lights do not turn on.	Power cable not securely connected. Remote control battery low or depleted. Power outlet malfunction.	Ensure all cables (mask to remote, remote to power adapter, power adapter to outlet) are firmly connected. Check if the remote control itself powers on. If not, replace its battery if applicable, or ensure it's receiving power. Test the power outlet with another device.
Mask feels uncomfortable.	Straps are too tight or improperly adjusted. Incorrect mask positioning.	Loosen and readjust the head straps for a comfortable fit. Ensure the mask is centered on your face, aligning eye and nose openings.
Lights flicker or dim.	Unstable power connection. Device malfunction.	Check all power connections. If the issue persists, contact customer support.

If you encounter issues not listed here or if solutions do not resolve the problem, please contact Pawfan customer support.

Specifications

Feature	Detail
Brand	Pawfan
Model Number	B0BWKCV2WL
UPC / GTIN	793288669648
Package Dimensions	16.26 x 9.96 x 2.05 inches
Weight	2.45 Pounds
Scent	Unscented
Skin Type	All
Special Feature	Lightweight, Portable

Warranty and Support

Specific warranty details are not provided in the product information. For warranty claims, technical support, or further assistance, please contact Pawfan customer service through their official website or the retailer where the product was purchased.

© 2025 Pawfan. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question