Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

AKIMRABY FTMT

AKIMRABY Bedside Commode Model FTMT User Manual

Brand: AKIMRABY | Model: FTMT

1. Introduction

This manual provides essential information for the safe and effective use of your AKIMRABY Bedside Commode Model FTMT. Designed for comfort and stability, this commode is suitable for elderly individuals, those with disabilities, and anyone requiring assistance with toileting. Please read all instructions carefully before assembly and use.



Image 1: AKIMRABY Bedside Commode Model FTMT in a bedroom setting.

2. Product Features

- **Robust & Safe Design:** Built for stability, supporting up to 500 lbs with a reinforced frame and non-slip base strips.
- **Detachable Armrests:** Offers a spacious 16.5-inch seat width and detachable armrests for adjustable comfort and easier transfers.
- **Adjustable Height:** Includes 4 additional extension legs to customize height from 16.5 to 18.9 inches.

- **Comfortable Use:** Features a soft PU seat cushion and a sponge backrest for sustained comfort and support.
- **Easy to Clean:** Comes with a 3.8-gallon inner bucket and is compatible with commode liners for simple waste disposal and minimal maintenance.

3. Assembly Instructions

Your AKIMRABY Bedside Commode is designed for straightforward assembly. Follow these steps to set up your commode:

1. **Base Assembly:** Ensure the main commode body is stable on a flat surface.
2. **Attach Height Extension Legs (Optional):** If additional height is desired, securely attach the four extension legs to the base of the commode. Ensure they click into place firmly.
3. **Install Detachable Armrests:** Align the armrests with the designated slots on the sides of the commode and push down until they lock securely.
4. **Place Backrest:** Position the sponge backrest onto the back of the commode.
5. **Insert Commode Bucket:** Place the 3.8-gallon inner bucket into the commode opening.
6. **Position Seat Cushion:** Place the soft PU seat cushion over the commode opening.



Image 2: Illustration of adjustable height options for personalized comfort.

4. Operating Instructions

4.1 Using Detachable Armrests

The commode features detachable armrests to facilitate easier transfers for users with limited mobility. To remove an armrest, locate the release mechanism (if present) or simply pull upwards firmly. To reattach, align and push down until it clicks into place.



Image 3: Detachable armrests provide a 16.5-inch seat width for user convenience.

4.2 Adjusting Height

The commode's height can be adjusted using the provided extension legs. The standard height is 16.5 inches, and with the extensions, it can be raised to 18.9 inches. Ensure all four legs are adjusted to the same height for stability.

4.3 Using the Commode Bucket and Liners

The commode includes a 3.8-gallon inner bucket for waste collection. For hygienic and easy disposal, it is recommended to use specially designed commode liners. Place a liner inside the bucket before use. After use, tie the liner securely and dispose of it according to local regulations.



Image 4: Using commode liners simplifies waste disposal and cleaning.

4.4 Comfort Features

The commode is equipped with a soft PU seat cushion and a sponge backrest to enhance user comfort, especially during extended periods of use.



Image 5: Soft PU seat cushion and sponge backrest for enhanced comfort.

4.5 Product Demonstration Video

Your browser does not support the video tag.

Video 1: A demonstration of the AKIMRABY Bedside Commode, showcasing its features and ease of use.

5. Cleaning and Maintenance

Regular cleaning ensures hygiene and extends the product's lifespan.

- **Commode Bucket:** After each use, remove the inner bucket, dispose of the liner (if used), and clean the bucket with a mild disinfectant and water. Rinse thoroughly and dry.
- **Commode Frame and Seat:** Wipe down the commode frame, seat, and backrest with a damp cloth and mild soap

or disinfectant. Avoid abrasive cleaners that may damage the material.

- **Non-slip Strips:** Periodically check the non-slip strips on the base for wear and tear to ensure continued stability.

6. Specifications

- **Product Dimensions:** 22 x 23 x 24 inches
- **Item Weight:** 13.85 Pounds
- **Model Number:** FTMT
- **Material:** PP, ABS, PU, PVC
- **Weight Capacity:** Up to 500 lbs
- **Seat Width:** 16.5 inches (with armrests)
- **Adjustable Height:** 16.5 to 18.9 inches
- **Bucket Capacity:** 3.8 gallons



Image 6: The commode is designed to support up to 500 lbs.

7. Safety Guidelines

- Always ensure the commode is placed on a flat, stable surface.
- Verify that all components, especially the height adjustment legs and armrests, are securely locked before each use.
- Do not exceed the maximum weight capacity of 500 lbs.
- Use caution when transferring to and from the commode. If necessary, seek assistance from a caregiver.
- Keep the commode clean to maintain hygiene and prevent the spread of germs.

8. Support & Warranty

For any questions, concerns, or warranty information regarding your AKIMRABY Bedside Commode Model FTMT, please contact the manufacturer directly. Refer to your product packaging or purchase documentation for specific contact details.

© 2024 AKIMRABY. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0CWCQNDQ6

More documents and user manuals: **https://manuals.plus/**