Klever Kits 60201

Klever Kits Microscope Explorer Set (Model 60201) - Instruction Manual

Introduction

Welcome to the world of microscopic discovery with your Klever Kits Microscope Explorer Set. This manual provides comprehensive instructions for assembling, operating, and maintaining your microscope, ensuring a safe and enriching experience. Designed for young scientists, this 2-in-1 microscope offers both stationary and portable observation modes, along with a phone holder for capturing your findings.

Package Contents

Before you begin, please verify that all components listed below are present in your package:



Figure 1: All components included in the Klever Kits Microscope Explorer Set.

- 1 Portable Microscope (60X-180X magnification)
- 1 Microscope Holder (base for stationary use)
- 1 Mobile Phone Holder
- 12 Prepared Specimen Slides (e.g., butterfly wing, dandelion fuzz)
- 4 Blank Slides
- 1 Lanyard
- 1 Cleaning Wipe
- 6 Anti-slip Sponge Mats
- 1 Light Source (requires 1 AAA battery, included for bottom light)

Setup

Follow these steps to set up your microscope for observation:

1. Battery Installation

The microscope requires one AAA battery for the bottom light source. Locate the battery compartment on the portable microscope unit and insert the battery according to the polarity markings. The top light is built-in and does not require separate battery installation.

2. Assembling the Microscope

The microscope can be used in two configurations: stationary with the holder or as a portable handheld unit.

Video 1: Demonstrates how to assemble the microscope, detach the portable unit, insert slides, and adjust focus.

1. **For Stationary Use:** Align the portable microscope unit with the microscope holder. Gently push the unit into the holder until it clicks securely into place. Ensure it is stable on a flat surface.
2. **For Portable Use:** To detach the portable microscope from its holder, press the body separation button located on the side of the holder and carefully pull the portable unit upwards.

3. Attaching the Phone Holder

The mobile phone holder allows you to capture images and videos through the microscope. Attach the phone holder to the eyepiece of the portable microscope unit. Secure your smartphone into the holder, ensuring the phone's camera lens is aligned with the microscope's eyepiece.

Operating Instructions

1. Using the Microscope (Stationary Mode)

1. **Place a Specimen Slide:** Gently slide a prepared specimen slide or a blank slide with your own sample onto the stage of the microscope holder. Ensure it is centered under the objective lens.
2. **Turn on the Light:** Use the switch to turn on the built-in LED light. You can choose between the upper light source (for opaque objects) or the bottom light source (for transparent slides).
3. **Adjust Magnification:** Rotate the zoom knob (labeled 60X-180X) to select your desired magnification level.
4. **Focus the Image:** Look through the eyepiece. Use the focus knob to bring the image into sharp view. If the image is still unclear, use the height adjustment knob (lever) on the side to move the slide stage up or down, then re-adjust the focus knob.



Figure 2: Labeled diagram of the stationary microscope components.

2. Using the Portable Microscope

Once detached from its holder, the portable microscope allows for on-the-go exploration.



Figure 3: A child using the portable microscope to examine a strawberry.

1. **Hold the Microscope:** Grip the portable microscope firmly.
2. **Position over Object:** Place the objective lens directly over the object you wish to observe (e.g., a leaf, fabric, skin).
3. **Adjust Magnification and Focus:** Use the zoom and focus knobs as described for stationary mode to achieve a clear image.



Figure 4: The portable microscope being used to observe a leaf outdoors.

3. Using the Phone Holder

To document your discoveries, use the mobile phone holder:



Figure 5: A smartphone attached to the microscope, showing a magnified view on its screen.

1. **Attach Phone:** Secure your smartphone in the phone holder, aligning its camera with the microscope's eyepiece.
2. **Open Camera App:** Launch your phone's camera application.
3. **Adjust View:** Adjust the microscope's focus and zoom until the magnified image is clear on your phone screen. You may need to slightly adjust the phone's position within the holder for optimal alignment.
4. **Capture:** Use your phone's camera controls to take photos or record videos of your observations.

Maintenance

- **Cleaning:** Use the provided cleaning wipe or a soft, lint-free cloth to gently clean the lenses and exterior of the microscope. Avoid abrasive cleaners or solvents.
- **Storage:** Store the microscope and its accessories in a clean, dry place away from direct sunlight and extreme temperatures. Keep slides in their protective case.
- **Battery Replacement:** Replace the AAA battery for the bottom light when it dims or stops working.

Troubleshooting

Problem	Possible Cause	Solution
Image is blurry or out of focus.	Incorrect focus or slide height.	Adjust the focus knob. If still blurry, use the height adjustment knob to move the slide stage up or down, then re-adjust the focus.
No light from the bottom source.	Dead battery or incorrect installation.	Check battery polarity or replace the AAA battery.
Cannot see anything through the eyepiece.	Lens cap on, slide not positioned correctly, or light off.	Ensure lens cap is removed. Center the slide. Turn on the light.
Phone camera not aligning with eyepiece.	Phone not properly secured or adjusted.	Carefully adjust the phone's position within the holder until the magnified view is centered on your screen.

Specifications

- **Model Number:** 60201
- **Magnification:** 60X - 180X
- **Light Source Type:** LED (Upper and Bottom)
- **Power Source:** 1 AAA Battery (for bottom light, included)
- **Material:** Plastic
- **Color:** White and Blue
- **Item Weight:** 14 ounces (0.4 Kilograms)
- **Package Dimensions:** 10.94 x 8.46 x 3.31 inches
- **Manufacturer:** JOVA US INC
- **UPC:** 840332675954

Safety Information

- This product is recommended for children aged 5 and above. Adult supervision is advised during use.
- Do not look directly into the LED light source as it may cause eye strain.
- Keep small parts, such as slides and batteries, away from young children to prevent choking hazards.
- Handle specimen slides with care to avoid breakage.
- Do not attempt to disassemble the microscope beyond the intended portable/stationary separation.

Warranty and Support

Klever Kits products are designed for quality and durability. For any questions, technical support, or warranty inquiries, please refer to the contact information provided on the product packaging or visit the official Klever Kits website. Please retain your proof of purchase for warranty claims.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question