Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

onerbl BP1001

onerbl AC-DC Adapter User Manual for Homiee BP1001 Digital Blood Pressure Monitor

Brand: onerbl | **Model Compatibility:** Homiee BP1001

1. Introduction

This document provides essential instructions for the safe and effective use of your onerbl AC-DC Adapter, designed as a replacement power supply for the Homiee BP1001 Digital Blood Pressure Monitor. This adapter converts standard AC wall power into the DC power required by your device, enabling operation or battery charging.

The adapter is manufactured with high-quality materials and has been tested to ensure compatibility and reliable performance with the specified device.



Image 1.1: The onerbl AC-DC Adapter with its integrated power cord. This image shows the compact design of the adapter unit and the neatly coiled cable.

2. Setup and Installation

Follow these steps to properly set up your AC-DC adapter:

1. **Verify Compatibility:** Ensure your device is the Homiee BP1001 Digital Blood Pressure Monitor or another device requiring the exact same power specifications.
2. **Connect to Device:** Insert the DC output plug of the adapter firmly into the power input port of your Homiee BP1001 Blood Pressure Monitor.
3. **Connect to Power Outlet:** Plug the AC input prongs of the adapter into a standard AC 100-240V wall power outlet.
4. **Power On:** Once connected, your device should receive power. Refer to your device's manual for specific power-on procedures.



Image 2.1: The onerbl AC-DC Adapter correctly plugged into a wall power outlet. This demonstrates the adapter's compact form factor when in use.

3. Operating Instructions

The onerbl AC-DC Adapter operates automatically once connected to a power source and your device. It provides a stable power supply for continuous operation or charging of your Homiee BP1001 Digital Blood Pressure Monitor.

- Ensure the adapter is securely plugged into both the wall outlet and your device.
- The adapter is designed for continuous use within its specified operating parameters.
- Disconnect the adapter from the wall outlet when not in use for extended periods.

4. Safety Information

To ensure safe operation and prevent damage, please observe the following safety guidelines:

- **Indoor Use Only:** This adapter is designed for indoor use in dry environments. Do not expose it to water, moisture, or high humidity.

- **Voltage Compatibility:** The adapter supports AC 100-240V, 50/60Hz input. Ensure your wall outlet voltage is within this range.
- **Overload Protection:** The adapter features built-in over voltage and over heat protection to safeguard both the adapter and your connected device.
- **Avoid Damage:** Do not bend, twist, or pull the cable excessively. Avoid placing heavy objects on the cable.
- **Ventilation:** Ensure adequate ventilation around the adapter during operation. Do not cover it with cloths or place it in enclosed spaces.
- **Disconnection:** Always grasp the adapter plug, not the cable, when disconnecting from the wall outlet.
- **Children and Pets:** Keep the adapter and cable out of reach of children and pets.



Image 4.1: Rear view of the onerbl AC-DC Adapter, displaying its model number, input voltage (100-240V AC, 50/60Hz), and the "CAUTION INDOOR USE ONLY" warning. The cable length is also indicated as 76 inches.

5. Maintenance and Care

Proper maintenance ensures the longevity and safe operation of your adapter:

- **Cleaning:** Disconnect the adapter from the power outlet and your device before cleaning. Use a soft, dry cloth to wipe the

surface. Do not use liquid cleaners or aerosols.

- **Storage:** When not in use, store the adapter in a cool, dry place, away from direct sunlight and extreme temperatures. Avoid tightly coiling the cable to prevent damage.
- **Inspection:** Periodically inspect the adapter and cable for any signs of damage, such as frayed wires, cracked casing, or bent prongs. If damage is found, discontinue use immediately.

6. Troubleshooting

If you encounter issues with your adapter, consider the following:

- **Device Not Powering On:**
  - Ensure the adapter is fully plugged into both the wall outlet and the device.
  - Test the wall outlet with another known working appliance to confirm it has power.
  - Verify that your device (Homiee BP1001) is functioning correctly with its original power source, if available.
- **Adapter Feels Hot:**
  - It is normal for adapters to feel warm during operation. However, if it is excessively hot to the touch, disconnect it immediately.
  - Ensure there is adequate ventilation around the adapter and it is not covered.
- **Visible Damage:** If the adapter or cable shows any signs of physical damage, discontinue use immediately and replace the unit.

7. Specifications

Product Type	AC-DC Power Adapter
Brand	onerbl
Compatibility	Homiee BP1001 Digital Blood Pressure Monitor
Input Voltage	AC 100-240V, 50/60Hz (Worldwide AC Input)
Output	Specific to Homiee BP1001 requirements (not explicitly stated in source data)
Protection Features	Over Voltage Protection, Over Heat Protection
Cable Length	Approximately 76 inches (193 cm)
Dimensions (Adapter Unit)	Approximately 3.1 inches (length) x 2.1 inches (width) x 1.2 inches (thickness)
Manufacturer	Shenzhen Waweis Technology Co., Ltd
ASIN	B0D9GHDGRK



Image 7.1: The onerbl AC-DC Adapter with its physical dimensions clearly marked: 3.1 inches in length, 2.1 inches in width, and 1.2 inches in thickness.

8. Warranty and Support

Specific warranty terms for this onerbl AC-DC Adapter are typically provided by the retailer or manufacturer at the time of purchase. Please refer to your purchase documentation for detailed warranty information.

For technical support or inquiries regarding this product, please contact the seller or manufacturer directly using the contact information provided with your purchase.

© 2024 onerbl. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0D9GHDGRK

More documents and user manuals: **https://manuals.plus/**