KingPavonini U8PRO

KingPavonini Wireless TENS Unit U8PRO Instruction Manual

Model: U8PRO | Brand: KingPavonini

1. Important Safety Information

Please read all instructions and warnings carefully before using this product. Failure to do so may result in injury or device malfunction.

- Consult a physician before use if you have any acute disease, tumor, infectious disease, heart disease, high fever, abnormal blood pressure, or skin allergies.
- Do not use if you are pregnant.
- Do not apply the device near the heart, on the head, face, throat, or any area with skin lesions or open wounds.
- Keep out of reach of children.
- Discontinue use immediately if you experience any discomfort or adverse reactions.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

2. Product Overview

The KingPavonini Wireless TENS Unit U8PRO is an innovative 3-in-1 muscle stimulator combining TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscle Stimulation), and massage capabilities for targeted pain relief and muscle relaxation. Its portable, dual-channel design allows for versatile use anywhere, anytime.

What's in the Box:

- 1 x Dual Channel Wireless TENS Machine (Charging Case)
- 2 x Wireless Main Units (Pucks)
- 4 x Large Rectangle Hydrogel Electrode Pads (3.23" x 1.86")
- 1 x USB-C Cable
- 1 x User Manual
- 1 x Instruction Card
- 1 x Gift Box

Key Features:

- **3-in-1 Functionality:** TENS, EMS, and Massage for comprehensive relief.
- **Wireless & Portable:** Compact design with a pocket-sized charging case for on-the-go use.

- **Dual Channel:** Allows simultaneous treatment of two different body areas.
- **6 Modes & 16 Intensity Levels:** Customizable settings for various massage techniques and pain relief needs.
- **Rechargeable Charging Case:** Stores and automatically recharges the main units.
- **Reusable Hydrogel Pads:** Flexible, comfortable, skin-friendly, and offer excellent self-adhesion for up to 45 uses.
- **Strong Magnetic Buckle:** Ensures a stable connection between the main unit and electrode pad.
- **LED Screen:** Clear display on the charging case for easy mode and intensity adjustment.



Image: KingPavonini Wireless TENS Unit U8PRO, showing the compact charging case and electrode pads.

3. Setup Guide

3.1 Charging the Device

1. Open the charging case. Place the two main units (pucks) into their respective slots inside the case.
2. Connect the provided USB-C cable to the charging port at the bottom of the case and to a suitable USB power source.
3. The LED screen on the charging case will indicate the charging status. A full charge typically takes 2 hours. The charging case holds a 450 mAh battery, providing extended use.



Image: The KingPavonini charging case connected via USB-C for recharging, illustrating the battery life of up to 12 hours with the case.

3.2 Preparing for Use

1. Ensure the main units are fully charged.
2. Take an electrode pad and remove the protective film from the adhesive side.
3. Align the main unit (puck) with the magnetic buckle on the electrode pad and attach it securely.
4. Clean and dry the skin area where you intend to apply the pad.
5. Apply the electrode pad with the attached main unit firmly to the desired skin area.



Image: Visual guide demonstrating the six steps for setting up and using the KingPavonini Wireless TENS Unit, from attaching the main unit to the pad to returning it to the charging case.

4. Operating Instructions

4.1 Power On/Off

- To turn on: Long press the power button on the main unit for 2 seconds. The LED indicator will light up.
- To turn off: Long press the power button on the main unit for 2 seconds. The LED indicator will turn off.

4.2 Mode Selection

- The device features 6 pre-set modes, including options like hammering and kneading.
- Press the 'M' button on the main unit or use the charging case's LED screen to cycle through the available modes.
- The LED screen on the charging case simplifies mode adjustment for clear operation.

4.3 Intensity Adjustment

- The device offers 16 intensity levels.
- Press the '+' button to increase intensity and the '-' button to decrease intensity on the main unit.
- Start with the lowest intensity and gradually increase until you feel a comfortable, yet effective, sensation.

4.4 Timer Function

- A 15-minute timer can be set for each session. The device will automatically shut off after the set time.

4.5 Dual Channel Operation

- The dual-channel design allows you to use both main units simultaneously on different body areas.
- Each main unit can be controlled independently for mode and intensity, offering personalized therapy.



Image: A woman using the KingPavonini Wireless TENS Unit on her back, highlighting the 6 manual modes and 16 intensity levels for personalized massage experience.

Your browser does not support the video tag.

Video: An official KingPavonini video demonstrating the wireless TENS unit, showcasing its compact design, ease of use, and various modes for muscle stimulation.

5. Maintenance and Care

5.1 Electrode Pad Care

- The hydrogel pads are reusable and washable. To maintain adhesion, gently clean the adhesive surface with a damp cloth or rinse with water after each use.
- Allow the pads to air dry completely before reattaching the protective film.
- Each pad can be reused up to 45 times with proper care. Replace pads when they lose their stickiness.



Image: A hand demonstrating the reusability of the hydrogel electrode pads, emphasizing their good stickiness and washable nature for up to 45 uses.

5.2 Device Storage

- Always store the main units and electrode pads in the charging case when not in use. This protects them from dust and damage and ensures they are charged for the next use.
- Store the device in a cool, dry place away from direct sunlight and extreme temperatures.

6. Troubleshooting

- **No Power/Weak Sensation:** Ensure the main units are fully charged. Check that the electrode pads are securely attached to the main units and firmly applied to clean, dry skin.
- **Pads Not Sticking:** Clean the adhesive surface of the pads with water and let them air dry. If adhesion does not improve, replace the pads.
- **Irregular Sensation:** Ensure both main units are properly connected and the pads are making full contact with the skin. Try adjusting the intensity level.

7. Specifications

Product Dimensions	2.2 x 1.82 x 3.7 inches
Item Weight	10.86 ounces
Model Number	U8PRO
Batteries	3 Lithium Polymer batteries (included)
Manufacturer	KingPavonini

8. Warranty and Support

KingPavonini offers a lifetime quality guarantee for this product. For any questions, concerns, or support needs, please contact our customer service team.
Customer Service: support@kingpavonini.com
Website: Visit the KingPavonini Store on Amazon

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question