Saved from manuals.plus.
Open the live page for links to other documents and user manuals.


PRO-SYS 5495-270

PRO-SYS JetFloss GO Cordless Water Flosser Instruction Manual

Model: 5495-270

Introduction

The PRO-SYS JetFloss GO is a cordless, rechargeable, and portable dental oral irrigator designed for effective plaque and debris removal. Its focused water jet refreshes and invigorates, contributing to improved oral hygiene and gum health. This manual provides essential information for the safe and proper use, maintenance, and care of your JetFloss GO device.



Effortless Flossing At Home or Away

LIGHTWEIGHT, COMPACT
Goes where you go, fits easily in your bag or luggage.

ADJUSTABLE WATER JET
Three settings for powering away debris and cleaning below the gumline.

5 QUICK CONNECT TIPS
Optimized for cleaning the hardest-to-reach places, including braces.

SUPERIOR 360° CONTROL
Swivels easily, so your hands twist less while flossing.



Image: The PRO-SYS JetFloss GO Cordless Water Flosser, a sleek black portable device.

Important Safety Information

Please read all instructions carefully before using this product.

- Keep away from children.
- Do not submerge the device in water.
- Charge the device only on a dry surface.
- If the device is damaged, do not attempt to repair it. Contact customer support for assistance.
- Dispose of the device and its components in accordance with local regulations.

Product Features

- **Cordless and Portable:** Designed for use at home or while traveling.
- **Rechargeable Battery:** Provides weeks of use on a single 4-hour charge.
- **IPX7 Waterproof:** Suitable for use in the shower.
- **Three Pressure Modes:** Adjustable water pressure from 40 to 100 PSI.
- **360° Rotating Tip:** Allows for easy cleaning of difficult-to-reach areas.

Setup and Initial Use

1. **Charge the Device:** Before first use, fully charge the JetFloss GO. Connect the charger cord to the device's charging port (ensure the port cover is open) and plug it into a power source. The battery charges in approximately four hours.
2. **Fill the Water Reservoir:** Open the reservoir cap and fill the reservoir with lukewarm water. Close the cap securely.
3. **Attach the Irrigator Tip:** Insert the desired irrigator tip into the opening on the top of the device until it clicks into place.

Operating Instructions

1. **Power On:** Press the power button to turn on the device.
2. **Select Pressure Mode:** Press the mode button to cycle through the three available pressure modes (Gentle, Standard, Pulse). Start with the lowest pressure setting and increase as comfortable.
3. **Position the Tip:** Lean over a sink and place the irrigator tip into your mouth, aiming it at your teeth and gumline.
4. **Begin Flossing:** Turn on the water jet. Guide the tip along your gumline, pausing briefly between teeth. Allow the water to flow into the sink.
5. **Complete Cleaning:** Continue until you have cleaned all areas of your mouth. The 360° rotating tip allows for easy access to all areas.
6. **Power Off:** Press the power button again to turn off the device.

Maintenance and Cleaning

- **After Each Use:** Empty any remaining water from the reservoir. Rinse the reservoir and irrigator tip under running water.

- **Deep Cleaning:** Periodically, clean the reservoir and tips with a mild soap solution. Ensure all soap residue is rinsed thoroughly.
- **Drying:** Allow all components to air dry before reassembling or storing the device.
- **Storage:** Store the JetFloss GO in a cool, dry place, away from direct sunlight. Use the provided travel bag for protection when traveling.

Troubleshooting

Problem	Possible Cause	Solution
Device does not turn on.	Low battery.	Charge the device fully.
Weak water pressure.	Low battery or clogged tip.	Charge the device or clean the irrigator tip. Ensure reservoir is full.
Water leaks from reservoir.	Reservoir cap not securely closed.	Ensure the reservoir cap is firmly closed.

If you experience issues not listed here, please contact PRO-SYS customer support.

Specifications

- **Model Number:** 5495-270
- **Product Dimensions:** 4.9 x 3.2 x 1.4 inches
- **Weight:** 13.4 ounces
- **Power Source:** Battery Powered (1 Lithium Ion battery included)
- **Water Pressure Range:** 40 - 100 PSI
- **Waterproof Rating:** IPX7
- **Manufacturer:** Pro-Sys

Warranty and Support

The PRO-SYS JetFloss GO comes with a **1-year limited warranty** from the date of purchase. This warranty covers manufacturing defects under normal use. It does not cover damage caused by misuse, accidents, or unauthorized modifications.
For warranty claims, technical support, or any questions regarding your device, please contact PRO-SYS customer service. Contact information can typically be found on the product packaging or the official PRO-SYS website.

Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

© 2026 PRO-SYS. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How do I get a replacement charging cable for my PRO-SYS JetFloss GO Cordless Water Flosser on the PRO-SYS 5495-270?

To obtain a replacement charging cable for your PRO-SYS JetFloss GO (Model 5495-270), you should contact the manufacturer directly as the manual does not list a specific part number for the cable alone. Here...

Read full answer

Open the live page

https://manuals.plus/asin/B0DJFRM3S7

More documents and user manuals: **https://manuals.plus/**