Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

MySmile LP233 Pro

MySmile LP233 Pro UVC Sterilizable Countertop Water Flosser Instruction Manual

Model: LP233 Pro

Introduction

The MySmile LP233 Pro UVC Sterilizable Countertop Water Flosser is designed to enhance your oral hygiene routine. This electric oral irrigator features a powerful pump and exclusive pulsation technology, delivering 1250-1650 pulses per minute and water stream strength of 50-180 PSI. It is clinically proven to remove up to 99.99% of plaque in treated areas and inhibit odor-causing bacterial growth.

Key features include multiple modes and pressure levels, a large water reservoir with a hydrogen-rich water filter, a UVC sterilizer light for jet tips, and a heating function for comfortable flossing.



Image: Visual representation of the water flosser's cleaning capability, showing removal of debris and plaque from teeth.

Important Safety Information

- Always unplug the device before cleaning.
- Do not use if you have a pacemaker or other electronic medical devices.
- Keep away from children.
- Do not look directly into the UV sterilizer light.

- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

Product Components

The MySmile LP233 Pro Water Flosser includes the following components:

- Main Unit (with 800ML water reservoir, control panel, and UVC sterilizer compartment)
- Magnetic Handle with pressure controls
- Variety of Jet Tips (Standard, Tongue Cleaner, Dental Plaque, Orthodontic, Periodontal)
- Hydrogen-rich Water Filter
- Power Cord




Image: An illustration showing the different types of jet tips included with the water flosser: Standard, Tongue Cleaner, Dental Plaque, Orthodontic, and Periodontal tips.

Setup Instructions

1. **Unpack the Device:** Carefully remove all components from the packaging.
2. **Place on a Stable Surface:** Position the main unit on a flat, stable countertop near an electrical outlet.
3. **Fill the Water Tank:** Open the water tank lid and fill the 800ml reservoir with lukewarm water. Ensure the hydrogen-rich water filter is properly installed.
4. **Attach a Jet Tip:** Select the desired jet tip and firmly insert it into the handle until it clicks into place.
5. **Connect Power:** Plug the power cord into the electrical outlet.



Image: A detailed view of the water flosser's 800ml water reservoir and the integrated hydrogen-rich water filter, highlighting its capacity and filtration system.

Operating Instructions

1. **Power On:** Press the ON/OFF button on the main unit to turn on the device.
2. **Select Mode:** Use the touch-sensitive panel to select one of the three modes:
    - **Child Mode:** 3 pressure levels, suitable for children (under adult supervision).
    - **Normal Mode:** 10 pressure levels, for daily and deep cleaning for adults.
    - **Pulse Mode:** 3 pressure levels, for gum massage and for individuals with mild periodontitis.
3. **Adjust Pressure:** Use the '+' and '-' buttons on the magnetic handle or the main unit to adjust the water pressure to your comfort level.
4. **Activate Heating Function (Optional):** Press the 'HEAT' button on the control panel if you desire warm water. The water will heat up to 20±2°C (68±4°F) in approximately 20 seconds.
5. **Begin Flossing:** Lean over the sink, place the jet tip in your mouth, and close your lips slightly to prevent splashing. Direct the water stream at a 90-degree angle to your gumline. Move the tip along your teeth and gumline, pausing briefly between teeth.
6. **Pause Function:** You can pause the water flow at any time by pressing the pause button on the magnetic handle.
7. **Power Off:** After flossing, press the ON/OFF button to turn off the device.




Image: The digital display and control panel of the water flosser, showing the three modes (Child, Normal, Pulse) and pressure level indicators.



Image: An internal view illustrating the heating element that warms the water for a more comfortable flossing experience.

Your browser does not support the video tag.

Video: An official product video demonstrating the features and usage of the MySmile LP233 Pro water flosser, including how to set it up and operate it for effective oral care.

Maintenance

1. **Empty and Dry Tank:** After each use, empty any remaining water from the reservoir and allow it to air dry to prevent mold and mildew.
2. **Clean the Main Unit:** Unplug the device. Wipe the exterior of the main unit with a soft, damp cloth. Do not immerse the unit in water.
3. **Sterilize Jet Tips:** After use, place the jet tips in the side storage compartment. The built-in UVC sterilizer light will automatically eliminate 99.99% of bacteria. The magnetic control ensures safety by shutting off power when the door is open.
4. **Replace Water Filter:** Periodically replace the hydrogen-rich water filter to ensure optimal performance and hygienic water. Refer to the filter packaging for replacement frequency.




Image: A visual of the UVC sterilizer light within the storage compartment, designed to eliminate bacteria from the jet tips.

Troubleshooting

If you encounter any issues with your MySmile LP233 Pro Water Flosser, please refer to the following common solutions:

- **No Water Flow:** Ensure the water tank is filled and the jet tip is securely attached. Check if the device is properly plugged in and powered on.
- **Weak Water Pressure:** Verify that the pressure level is set appropriately for your needs. Clean the jet tip to ensure no blockages.
- **Device Not Turning On:** Check the power connection. Ensure the power outlet is functional.
- **UVC Sterilizer Not Working:** Ensure the storage compartment door is fully closed. The magnetic control will prevent operation if the door is open.

For further assistance, please contact MySmile customer support.

Specifications

Model Number	LP233 Pro
Brand Name	MySmile
Tank Volume	800 Milliliters (27 Ounces)
Water Pressure Range	50-180 PSI
Pulsation Rate	1250-1650 pulses/minute
Operation Modes	Child (3 levels), Normal (10 levels), Pulse (3 levels)
Heating Function	Heats water to 20±2°C (68±4°F) in approx. 20 seconds
Sterilization	UVC Sterilizer Light (eliminates 99.99% bacteria)
Power Source	Corded Electric
Material Type	Plastic, Metal, Glass, Silicone
Certifications	CE, FCC, RoHS
Package Dimensions	9.61 x 8.94 x 7.2 inches
Item Weight	4.54 Pounds

Warranty and Support

For information regarding product warranty, returns, or technical support, please contact MySmile customer service. You can typically find contact details on the product packaging or the official MySmile website.
Manufacturer: MySmile

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0DJX3K95H

More documents and user manuals: **https://manuals.plus/**