Ivifloae Ivifloaesafxv97tub

Ivifloae Pill Counting Tray with Spatula Instruction Manual

Model: Ivifloaesafxv97tub

Introduction

Thank you for choosing the Ivifloae Pill Counting Tray with Spatula. This manual provides essential information for the safe and effective use of your new tablet counting device. This product is designed to facilitate the hygienic and efficient counting and dispensing of tablets and capsules in various settings, including home and professional pharmacy environments. Please read these instructions carefully before first use and retain them for future reference.

Product Overview and Components

The Ivifloae Pill Counting Tray consists of a main counting tray and a small spatula. It is constructed from durable, odor-free plastic, ensuring a hygienic and long-lasting tool for medication management.

- **Counting Tray:** The main blue plastic surface with a clear hinged cover. It features a lower left channel for dispensing and an upper right channel for returning excess tablets.
- **Spatula:** A small, flat tool used to manipulate and count tablets on the tray without direct hand contact.
- **Hinged Cover:** A transparent cover that protects the tablets and facilitates pouring.



Image 1: Angled view of the Ivifloae Pill Counting Tray with its clear hinged cover and integrated spatula.



Image 2: Close-up view of the durable hinge mechanism connecting the clear cover to the main tray.

Setup

The Ivifloae Pill Counting Tray requires minimal setup. Before first use, and after each use, ensure the tray and spatula are clean and dry.

1. Unpack the counting tray and spatula from its packaging.
2. Wash the tray and spatula with mild soap and warm water.
3. Rinse thoroughly and allow to air dry completely before use.
4. Place the tray on a clean, flat, and stable surface.

Operating Instructions

Follow these steps to accurately count and dispense tablets using your Ivifloae Pill Counting Tray:

1. **Prepare the Tray:** Ensure the tray is clean and dry. Place it on a stable surface.
2. **Dispense Tablets onto Tray:** Carefully pour the required amount of tablets or capsules from their container onto the main

counting surface of the tray. Avoid overfilling.

3. **Count Tablets:** Use the provided spatula to gently spread the tablets across the tray. Systematically move and count the tablets, using the spatula to guide them into groups or lines as needed. The design allows for counting without direct hand contact, promoting hygiene.
4. **Dispense into Container:** Once the desired number of tablets has been counted, guide them towards the lower left channel of the tray. Tilt the tray slightly to allow the tablets to slide smoothly through the channel into your target container (e.g., a daily pill box or prescription bottle).
5. **Return Excess Tablets:** If there are any remaining tablets on the tray after dispensing, guide them towards the upper right channel. Tilt the tray to pour these excess tablets back into their original storage bottle.



Image 3: Illustration of pills being poured onto the counting tray, ready for counting with the spatula.



Image 4: The tray being used to dispense counted pills into a medication bottle via the lower left channel.

Maintenance

Proper maintenance ensures the longevity and hygienic performance of your counting tray.

- **Cleaning:** After each use, wash the tray and spatula with mild soap and warm water. Avoid abrasive cleaners or harsh chemicals that could damage the plastic.
- **Drying:** Rinse thoroughly and allow all components to air dry completely before storing. Ensure no moisture remains, especially in crevices.
- **Storage:** Store the clean and dry tray and spatula in a clean, dry place, away from direct sunlight and extreme temperatures.

Troubleshooting

Problem	Possible Cause	Solution

Problem	Possible Cause	Solution
Tablets stick to the tray.	Static electricity or residue from previous use.	Ensure the tray is completely dry and clean. A very light wipe with an anti-static cloth (if safe for plastic) might help, but usually, proper cleaning is sufficient.
Difficulty guiding tablets with spatula.	Incorrect angle or too many tablets on the tray.	Use a gentle, sweeping motion. Do not overload the tray; count in smaller batches if necessary. Practice will improve technique.

Specifications

Model Number	lvifloaesafxv97tub
Material	Plastic
Dimensions (Approx.)	8.27 x 6.69 x 1.57 inches (21 x 17 x 4 cm)
Weight (Approx.)	5.86 ounces (166 grams)
Color	Blue (Tray), Clear (Cover)

Safety Information

- Keep out of reach of children and pets.
- This product is intended for counting and dispensing tablets and capsules only. Do not use for other purposes.
- Always ensure the tray is clean and dry before and after each use to prevent cross-contamination.
- Handle with care to avoid dropping or damaging the tray, which could affect its functionality.

Warranty and Support

For information regarding product warranty or customer support, please refer to the retailer's return policy or contact the seller directly. Please retain your proof of purchase for any warranty claims.
Seller: beiyuanwangluo
Amazon Link: Seller Information

© 2024 lvifloae. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question