Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

Oricle HA-JHA490-2.0

Oricle Standard 2.0 Hearing Aids User Manual

Model: HA-JHA490-2.0 | Brand: Oricle

1. Introduction

Welcome to the Oricle Standard 2.0 Hearing Aids user manual. These advanced hearing aids are designed to provide clear sound and comfortable wear for adults, featuring advanced noise cancellation and easy volume control. This manual will guide you through the setup, operation, and maintenance of your new hearing aids to ensure optimal performance and longevity.



Image: Oricle Standard 2.0 Hearing Aids with their portable charging case, showcasing the compact design.

2. What's in the Box

Upon opening your Oricle Standard 2.0 Hearing Aids package, you should find the following items:

- Oricle Standard 2.0 Hearing Aids (Left and Right)
- Portable Charging Case
- Charging Cable (Charger)
- Multiple Sizes of Ear Tips (for comfortable fit)
- User Manual (this document)

3. Product Features

The Oricle Standard 2.0 Hearing Aids are equipped with several key features designed to enhance your hearing experience:

- **No-More-Ear-Discomfort:** Experience all-day comfort with these in-ear silicone earplugs. Enjoy a secure, gentle fit with soft ear tips for extended wear.
- **Clear Conversations in Any Environment:** This OTC hearing aid enhances speech clarity, making it easier to hear family and friends even in noisy settings. Enjoy confident, effortless connections in every environment.
- **Premium-Comfort-Design:** This completely-in-canal hearing aid for adults features a sleek, compact design that remains stable and comfortable for long-term wear, without concerns about slipping.
- **Exceptional Battery Life:** These CIC hearing aids work up to 60 hours on a single charge. The portable charging case lets you conveniently recharge your hearing aids anytime, whether you're at work, traveling, or on the move.
- **Easy-to-Use Controls:** Long press for 5 seconds to turn On/Off. Features Auto on/off when taken out or placed into the compact charging case with smart memory technology to keep initial and desired settings. Choose from 5 different sound levels for adapting environments.



Image: Visual representation of the Oricle hearing aids highlighting their invisible design, clear sound, digital core, long battery life, and compact case.



Image: Depiction of the versatility of Oricle hearing aids in different environments such as concerts, nature, planes, watching TV, phone calls, and restaurants.



Image: Infographic detailing the sound quality features including Hi-Fi, Ultra Crystal Clear, HD Surround, Intelligent Noise Reduction, and Zero Whistling.

4. Setup Guide

4.1 Initial Charging

Before first use, ensure your Oricle hearing aids and charging case are fully charged. Place the hearing aids into their respective slots in the charging case. Connect the charging cable to the case and a power source. The indicator lights on the case will show charging status. A full charge takes approximately 2 hours.



Image: Details on the premium rechargeable battery, showing 2 hours for a full charge, 16 hours of working time, and 80 hours of total power with the charging case.

4.2 Selecting Ear Tips

The Oricle Standard 2.0 Hearing Aids come with multiple sizes of silicone ear tips. Select the size that provides a secure, comfortable fit in your ear canal. A proper fit is crucial for optimal sound quality and to prevent feedback.

4.3 Insertion and Removal

Gently insert the hearing aid into your ear canal. The discreet design allows for a comfortable and virtually unnoticeable fit. To remove, gently pull on the small tab extending from the hearing aid.



Image: A visual comparison demonstrating the comfortable and discreet fit of Oricle hearing aids, suitable for wearing with glasses and masks, unlike some other designs.

5. Operating Instructions

5.1 Power On/Off

- **Automatic:** The hearing aids will automatically power on when removed from the charging case and power off when placed back into the case.
- **Manual:** If not using the charging case, long press the button on the hearing aid for 5 seconds to manually turn it On or Off.

5.2 Volume Control

The Oricle Standard 2.0 Hearing Aids feature 5 different sound levels to adapt to various environments. To adjust the volume, press the button on the hearing aid. Each press will cycle through the available sound levels.

5.3 Smart Memory Technology

The hearing aids are equipped with smart memory technology, which retains your initial and desired settings. This means you won't need to re-adjust settings every time you use them after charging.



Image: Diagram illustrating the easy-to-use design, including auto on/off with the charging case, smart memory, manual on/off via a 5-second long press, and 5 sound levels.

5.4 Official Product Video

Your browser does not support the video tag.

Video: An official product video demonstrating the Oricle 2.0 Hearing Aids, including their design, features, and ease of use. (Duration: 0:36)

6. Maintenance

Proper maintenance will extend the life of your Oricle hearing aids:

- **Cleaning:** Regularly clean the ear tips and the body of the hearing aids with a soft, dry cloth. Do not use water or cleaning solutions directly on the electronic components.
- **Storage:** When not in use, store the hearing aids in their charging case to protect them from dust and damage, and to keep them charged.

- **Battery Care:** Avoid completely draining the battery frequently. Recharge as needed. The portable charging case provides convenient recharges on the go.

7. Troubleshooting

If you encounter any issues with your Oricle Standard 2.0 Hearing Aids, please refer to the following common troubleshooting tips:

Problem	Possible Cause	Solution
No Sound / Low Sound	Low battery; Ear tip clogged with earwax; Incorrect volume setting.	Charge the hearing aids; Clean the ear tip; Adjust volume level.
Charging Case Light Stays Red	Case not fully charged; Connection issue.	Ensure the charging cable is securely connected to both the case and a working power source. Allow sufficient time for a full charge (approx. 2 hours).
Whistling or Feedback Noise	Improper fit; Earwax buildup in ear canal.	Try a different size ear tip for a snugger fit; Consult a healthcare professional for earwax removal if necessary.
Hearing Aid Not Turning On/Off Automatically	Improper placement in charging case; Sensor obstruction.	Ensure hearing aids are correctly seated in the charging case slots. Clean any debris from the charging contacts.

For more detailed troubleshooting tips, please refer to the comprehensive manual included in your product package or contact Oricle customer support.

8. Specifications

Feature	Detail
Model Name	Standard 2.0 Hearing Aids
Item Model Number	HA-JHA490-2.0
Brand	Oricle
Fit Type	In-Ear (CIC - Completely-in-Canal)
Special Features	Lightweight, Noise Reduction, Rechargeable, Volume Control
Battery Type	Lithium Ion (1 included)
Battery Weight	43 Grams
Approx. Battery Life (Single Charge)	Up to 60 hours
Charging Case Total Power	Up to 80 hours (multiple recharges)
Package Dimensions	4.92 x 4.02 x 2.2 inches

Feature	Detail
Item Weight	7.83 ounces
Manufacturer	Oricle
First Available Date	November 7, 2024

9. Warranty and Support

For information regarding product warranty, returns, or technical support, please refer to the warranty card included in your product packaging or visit the official Oricle website. You may also contact Oricle customer service directly for assistance.

Manufacturer: Oricle

© 2025 Oricle. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0DMBX8RZT

More documents and user manuals: **https://manuals.plus/**