Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

AKIMRABY B0DNQR5MTR

AKIMRABY Heavy Duty Bedside Commode Instruction Manual

Model: B0DNQR5MTR

1. Product Overview

The AKIMRABY Heavy Duty Bedside Commode is a versatile and robust mobility aid designed to provide comfort and convenience for elderly individuals, disabled persons, postoperative patients, and pregnant women. This commode chair features a sturdy design with a high weight capacity, adjustable height options, and detachable armrests for ease of use and transfer. It seamlessly transitions between a bedside commode and a comfortable sofa chair, making it a practical solution for various needs.
The complete care package includes 120 commode liners and 50 large commode pads for easy waste disposal and maintenance.

Video: Overview of the AKIMRABY Bedside Commode, showcasing its dual functionality as a commode and a soft cushion chair, its robust 900 lbs load capacity, and features like reinforced armrests and adjustable height for user support and comfort.

2. Setup and Assembly

Follow these steps to assemble your AKIMRABY Bedside Commode. For a visual guide, please refer to the installation video below.

1. **Attach the Legs:** Place the main body of the commode upside down. Insert the four main legs into their corresponding slots. Secure them firmly using the provided screws.
2. **Install Extension Legs (Optional for Height Adjustment):** If additional height is desired, install the extension legs onto the main legs. The commode offers three height options: 11.5" (low), 15.7" (standard), and 17.7" (elevated).
3. **Attach Anti-Slip Strips:** For enhanced stability and floor protection, attach the anti-slip strips to the bottom of each leg.
4. **Insert the Backrest:** Position the backrest into the designated slot on the main body. Secure it with the provided fixed dowels.
5. **Attach Detachable Armrests:** Attach the armrests to the sides of the main body. These can be easily removed to facilitate transfers for users with limited mobility.
6. **Install Hanging Tissue Box:** The hanging tissue box can be installed on either side of the commode for convenient access to tissues or other small items.
7. **Place Inner Bucket and Cushion:** Insert the 3.6-gallon inner commode bucket into the main body. Place the soft cushion over the toilet seat for comfortable seating when used as a chair.



Image: The AKIMRABY Bedside Commode illustrating its dual functionality: a bedside commode with a toilet seat and a soft cushion chair. This highlights its versatile design for various user needs.



Image: A detailed view of the commode's sturdy structure, emphasizing its 4-large support legs that provide a 900 lbs load capacity. It also shows the 17-inch seat width for comfortable seating.



Image: Illustration of the detachable armrests, showing how they can be removed to assist users with limited mobility in transferring to and from the commode.



Image: Visual representation of the three adjustable height options: 11.5" (low), 15.7" (standard), and 17.7" (elevated), achieved by adding or removing extension legs.

Video: Detailed installation guide for the AKIMRABY Bedside Commode, demonstrating step-by-step assembly of the legs, backrest, and armrests.

3. Operating Instructions

3.1 Using as a Bedside Commode

- **Positioning:** Place the commode next to the bed or in a convenient location. Ensure the commode is stable on a flat surface.
- **Height Adjustment:** Adjust the commode height to match the user's needs for comfortable and safe transfer. Use the extension legs to achieve 11.5", 15.7", or 17.7" height.
- **Detachable Armrests:** For easier transfer from a bed or wheelchair, the armrests can be detached by pressing the release buttons and lifting them off. Reattach them securely after transfer for added support.
- **Commode Bucket and Liners:** Lift the toilet seat and the inner cover to access the 3.6-gallon commode bucket. Place a commode liner inside the bucket, ensuring it covers the edges. The included 120 liners and 50 large commode pads are designed for easy waste disposal.
- **Slow-Close Design:** The toilet seat features a slow-close mechanism to prevent sudden slamming, ensuring quiet operation, especially during night use.
- **Sealing:** The commode bucket includes a cover that provides excellent sealing to contain odors, making it suitable for indoor use.

3.2 Using as a Sofa Chair

- **Cover the Commode:** Ensure the commode bucket is covered, and the toilet seat is lowered.
- **Place Cushion:** Place the soft cushion over the toilet seat to transform the commode into a comfortable chair.
- **Backrest Support:** The ergonomic backrest provides strong support for comfortable seating.

4. Maintenance and Cleaning

- **Waste Disposal:** After use, carefully remove the commode liner from the bucket, tie it securely, and dispose of it according to local regulations. The included commode pads convert liquid waste into a gel for easier and cleaner disposal.
- **Cleaning the Bucket:** Regularly clean the 3.6-gallon commode bucket with a mild disinfectant and water. Rinse thoroughly and allow it to dry before placing a new liner.
- **Cleaning the Chair:** Wipe down the commode chair's surfaces, including the seat, backrest, and armrests, with a damp cloth and mild soap. Avoid abrasive cleaners that may damage the material.
- **Anti-Slip Strips:** Periodically check the anti-slip strips on the legs to ensure they are clean and intact for optimal stability.



Image: Close-up of the enlarged commode pads, illustrating their super absorption feature that converts 70 oz of liquid into a gel within seconds for

convenient and hygienic disposal.

5. Troubleshooting

If you encounter any issues with your AKIMRABY Bedside Commode, please review the assembly instructions to ensure all components are correctly installed. For persistent problems or concerns, please contact AKIMRABY customer support for assistance.

6. Specifications

- **Model:** B0DNQR5MTR
- **Weight Capacity:** Up to 900 lbs
- **Seat Width:** 17 inches
- **Height Options:** 11.5" (low), 15.7" (standard), 17.7" (elevated)
- **Commode Bucket Capacity:** 3.6 gallons
- **Materials:** Reinforced legs, anti-slip strips, soft cushion, durable plastic construction
- **Included Accessories:** 120 commode liners, 50 large commode pads, detachable armrests, hanging tissue box

7. Warranty and Support

Specific warranty information for the AKIMRABY Bedside Commode (Model B0DNQR5MTR) is not provided within this manual. For detailed warranty terms, conditions, and to register your product, please refer to the product packaging or contact AKIMRABY customer support directly. For any support inquiries, technical assistance, or replacement parts, please visit the official AKIMRABY website or use the contact information provided on your product's packaging.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0DNQR5MTR

More documents and user manuals: **https://manuals.plus/**