Serianie 0005

Serianie Model 0005 Electric Wheelchair User Manual

Your Guide to Safe and Efficient Operation

Introduction Safety Precautions Setup & Assembly Operating Guidelines Maintenance & Care Specifications Support

1. Introduction

This user manual provides essential information for the safe and effective operation, maintenance, and care of your Serianie Model 0005 Electric Wheelchair. Please read this manual thoroughly before using the wheelchair to ensure proper understanding of its features and functions. Keep this manual for future reference.



Figure 1: Serianie Model 0005 Electric Wheelchair

2. Safety Precautions

Your safety is paramount. Adhere to all safety guidelines to prevent injury or damage to the wheelchair.

- Always ensure the wheelchair is fully charged before extended use.
- Familiarize yourself with the joystick controls and speed settings in a safe, open area before operating in crowded or complex environments.
- Use the safety belt provided to secure yourself in the chair during operation.
- Avoid operating the wheelchair on excessively steep inclines or unstable terrain beyond its specified capabilities.
- Regularly inspect tires, brakes, and electrical connections for any signs of wear or damage.
- Do not attempt to modify the wheelchair or its components, as this may compromise safety and void the warranty.

3. Package Contents

The package includes the following items:

- One lightweight folding electric wheelchair designed for adults.
- Removable battery pack.
- Joystick controller.
- Charging cable and adapter.
- User manual (this document).

4. Setup and Assembly

The Serianie Model 0005 Electric Wheelchair is designed for easy setup. Follow these steps for initial assembly and preparation:

1. **Unfolding the Wheelchair:** Carefully remove the wheelchair from its packaging. Place it on a flat surface and gently unfold the frame until it locks into place. Ensure all locking mechanisms are secure.
2. **Battery Installation:** Locate the battery compartment, typically beneath the seat. Insert the removable battery pack, ensuring it clicks securely into position.
3. **Connecting the Joystick:** Attach the joystick controller to the designated port on the armrest. Ensure the connection is firm.
4. **Footrest Adjustment:** Adjust the footrests to a comfortable height for the user.



Figure 2: Folding the Serianie Electric Wheelchair



Figure 3: Compact folded size for storage

Your browser does not support the video tag.

Video 1: Demonstrating the folding and unfolding of the electric wheelchair.

5. Operating Guidelines

The Serianie Model 0005 is designed for intuitive control and versatile use, both indoors and outdoors.

5.1. Joystick Control

The wheelchair is controlled via a 360° joystick. Push the joystick in the desired direction to move. The further you push, the faster the wheelchair will move, up to its maximum speed.



Figure 4: Joystick Controller with Speed Settings

5.2. Speed Adjustment

The wheelchair features 5 adjustable speed levels. Use the speed adjustment buttons on the joystick to increase or decrease your speed. Start with a lower speed setting until you are comfortable with the controls.



Figure 5: 5 Levels Adjustable Speed

5.3. All-Terrain Capability

This electric wheelchair is designed for multi-terrain application, providing a smooth ride on various surfaces. However, exercise caution on uneven or soft ground.



Figure 6: All-Terrain Use



Figure 7: Navigating various outdoor environments

6. Maintenance and Battery Care

Proper maintenance ensures the longevity and optimal performance of your electric wheelchair.

- **Battery Charging:** Connect the charging cable to the battery port and plug the adapter into a standard electrical outlet. The battery takes approximately 5 hours for a full charge.
- **Battery Capacity:** The battery has a capacity of 12 Ah, providing an average travel range of 25 miles on a full charge.
- **Cleaning:** Wipe down the wheelchair frame and seat with a damp cloth. Avoid using harsh chemicals or abrasive cleaners.
- **Storage:** When not in use, store the wheelchair in a dry, cool place. If storing for an extended period, ensure the battery is fully charged and periodically recharge it to maintain battery health.

Figure 8: Upgraded Removable Battery Details

7. Troubleshooting

If you encounter any issues with your Serianie Model 0005 Electric Wheelchair, please refer to the following common troubleshooting steps:

- **Wheelchair not powering on:** Check if the battery is properly installed and charged. Ensure the power button on the joystick is pressed.
- **Wheelchair not moving:** Verify that the joystick is securely connected. Check for any obstructions around the wheels. Ensure the brakes are disengaged (if applicable).
- **Unusual noises:** Stop using the wheelchair immediately and inspect for any loose parts or foreign objects caught in the wheels or motor. If the noise persists, contact customer support.

For issues not covered here, please contact Serianie customer support.

8. Specifications

Feature	Detail
Product Dimensions	41.3 x 24.8 x 37 inches (D x W x H)
Item Weight	100 Pounds
Maximum Weight Recommendation	300 Pounds
Seat Width	18.5 Inches
Seat Height	21.6 Inches
Seat Depth	17.3 inches
Wheel Size	16 Inches
Turn Radius	4 Feet
Brake Style	Electromagnetic
Maximum Speed	6 Kilometers per Hour
Power Source	Battery Powered
Battery Average Life	25 Miles
Wheel Type	Pneumatic
Frame Material	Alloy Steel
Color	Black
Model Number	0005



Figure 9: Seat Dimensions




Figure 10: Payload Capacity

9. Warranty and Customer Support

Serianie stands by the quality of its products. Your electric wheelchair comes with a warranty. For any inquiries, technical support, or issues regarding your electric wheelchair, please contact our customer service team. We are committed to helping you solve any problem within 24 hours.
For further assistance, please refer to the contact information provided with your purchase documentation.

© 2026 Serianie. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question