Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

SALLOOT DT015 Pro

SALLOOT IPL Laser Hair Removal Device

MODEL: DT015 PRO - USER MANUAL

1. Introduction

Thank you for choosing the SALLOOT IPL Laser Hair Removal Device, Model DT015 Pro. This advanced device is designed to provide effective and comfortable at-home hair removal using Intense Pulsed Light (IPL) technology. Featuring dual lights, sapphire ice cooling, and a powerful 4-pulse PHR mode, the DT015 Pro offers a salon-grade experience for both women and men. Please read this manual thoroughly before first use to ensure safe and optimal operation.

2. Safety Information

Your safety is paramount. Adhere to the following guidelines:

- Always wear the provided safety goggles during operation to protect your eyes from the intense light pulses.
- Perform a patch test on a small skin area 24 hours before full treatment to check for adverse reactions.
- Ensure your skin tone and hair color are compatible with IPL technology. This device is most effective on light to medium skin tones with dark hair. It is not suitable for very dark skin tones or very light (blonde, red, white, gray) hair colors. Refer to the skin tone chart in the 'Operating Instructions' section.
- Do not use on tattooed areas, moles, birthmarks, or open wounds.
- Keep the device dry and do not use it in wet environments.
- Keep out of reach of children.

For any concerns regarding applicability or usage, please contact SALLOOT customer support.

Ingredients: Sapphire, Quartz lamp, Glass, Resin

3. Product Components

The SALLOOT IPL Laser Hair Removal Device package includes the following items:

- SALLOOT IPL Laser Hair Removal Device (DT015 Pro)
- Power Adapter and Cord
- Safety Goggles
- Razor
- Cleaning Cloth
- User Manual



Figure 3.1: SALLOOT IPL Laser Hair Removal Device with accessories.

4. Setup Guide

Follow these steps to prepare your device for use:

1. **Prepare Skin:** Shave the desired treatment area thoroughly. Ensure the skin is clean, dry, and free from any lotions, creams, or deodorants.
2. **Connect Power:** Plug the power adapter into the device and then into a power outlet.
3. **Wear Safety Goggles:** Put on the provided safety goggles to protect your eyes from the intense light flashes.
4. **Download AI Coaching App (Optional but Recommended):** The i15 Pro features an AI Coaching APP that auto-recommends energy levels and flash modes based on your skin tone, hair density, and treatment area. Scan the QR code in the user manual or search for "CTONE" in your app store (Google Play or App Store) to download. The app provides personalized treatment plans, step-by-step tutorials, and smart email reminders.




Figure 4.1: Utilizing the AI Coaching App for personalized treatment.

5. Operating Instructions

Understand how to effectively use your SALLOOT IPL device:

5.1 Skin Tone and Hair Color Compatibility

IPL technology targets melanin in hair follicles. For effective and safe treatment, ensure your skin tone and hair color are compatible:




Figure 5.1: Skin Tone and Hair Color Compatibility.

5.2 Power On/Off

Press and hold the power button to turn the device on or off.

5.3 Mode Selection and Features

- **Dual Lights 100W:** The Trinity i15 Pro utilizes an industry-leading 100W high power and a 4.2cm² light window with

dual lights, ensuring 4X efficiency and fast results.

- **4-Pulse Burst PHR Mode:** This mode boasts a total energy output of 27J per 4-pulse burst, designed to eliminate stubborn and ingrown hair by carbonizing the hair follicle up to 98%.
- **Sapphire Ice Cooling:** The enlarged 16.6cm² sapphire ice-cooling window instantly cools the contact surface to 35°F (2°C) and maintains a 53°F (11.7°C) ice-cooling touch for up to 60 minutes of continuous use, ensuring comfort, especially for sensitive areas.
- **Auto Gliding Fast Mode:** For large areas like legs and arms, long press the flashing button for 2 seconds to switch to auto gliding mode for continuous flashes.
- **Manual Flash Mode:** For smaller or sensitive areas such as the face, bikini line, or underarms, use the manual flash mode for precise control.



Figure 5.2: PHR Mode and Sapphire Ice Cooling.

5.4 Energy Levels

Choose from 15 sensitivity levels to adjust the intensity according to your skin's sensitivity and the specific area being treated. Start with a lower level and gradually increase as your skin tolerates.

5.5 Treatment Areas

The device is suitable for whole-body treatment, including legs, arms, bikini line, face, armpit, back, chest, fingers, and toes.



Figure 5.3: Whole Body Treatment Areas.

5.6 Usage Frequency

Use the device multiple times based on your hair growth. Pairing it with the smart app is recommended for optimal results.

- **Weeks 1-2:** 3-4 times/week
- **Weeks 3-4:** 2-3 times/week
- **After 4 weeks:** 2 times/monthly for long-lasting results




Figure 5.4: Recommended Usage Frequency.

6. Maintenance

Proper maintenance ensures the longevity and effectiveness of your device:

- **Cleaning:** After each use, gently wipe the light window and the device surface with the provided cleaning cloth to remove any hair or residue. Do not use abrasive cleaners or harsh chemicals.
- **Storage:** Store the device in its original packaging or a protective case in a cool, dry place, away from direct sunlight and moisture.

7. Troubleshooting

If you encounter any issues with your SALLOOT IPL device, refer to the following common solutions:

- **Device Overheating/Automatic Shutdown:** If the device overheats and automatically shuts down, particularly when using the PHR quad-pulse mode continuously for more than 10 minutes, unplug the device, wait 10 seconds, and then plug it back in. Allow it to cool down before resuming use.
- **Device Not Turning On:** Ensure the power adapter is securely connected to both the device and a working power outlet.
- **Ineffective Hair Removal Results:**
  - Verify your skin tone and hair color are compatible with IPL technology (refer to Section 5.1).
  - Ensure the light window is in full contact with your skin during flashes.
  - Increase the energy level if comfortable and appropriate for your skin type.
  - Maintain consistent usage according to the recommended frequency (refer to Section 5.6).

If the issue persists, please contact customer support.

8. Specifications

Feature	Detail
Model Number	DT015 Pro
Package Dimensions	10.31 x 6.54 x 4.69 inches
Item Weight	3.62 Pounds
UPC	642380741298
Manufacturer	SALLOOT
Optical Outlet	4.2cm²
Ice-cooling Area	16.6cm²
Flash Speed	0.15s per flash
Flashes per Minute	400 flashes/minute

9. Warranty and Support

SALLOOT is committed to providing excellent customer service and product reliability.

- **Extended Warranty:** Your purchase may include an extended warranty. Please refer to the product registration information or contact customer support for details.
- **Lifetime 24/7 Support:** SALLOOT offers lifetime 24/7 support for its products. You can join their community or contact them directly for assistance.
- **Contact Information:**
  Email: support@sallootbeauty.com
  Website: www.sallootbeauty.com

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

How can I resolve the issue where the SALLOOT DT015 Pro flashlight fails to trigger in both manual and automatic modes?

If your SALLOOT DT015 Pro is powered on and the cooling fan is running, but the device is not flashing in either Manual or Auto Gliding mode, please check the following troubleshooting steps: Skin...

Read full answer

Open the live page

https://manuals.plus/asin/B0DPM7XXR7

More documents and user manuals: **https://manuals.plus/**