Jaioyios AF1378-1

Jaioyios Electric Chair Lift Instruction Manual

Model: AF1378-1
Brand: Jaioyios

1. Introduction and Overview

This manual provides essential instructions for the safe and effective use of your Jaioyios Electric Chair Lift, Model AF1378-1. Please read this manual thoroughly before assembly, operation, or maintenance to ensure proper function and user safety. This device is designed to assist individuals in rising from the floor and facilitating transfers within an indoor environment.



The Jaioyios Electric Chair Lift is designed for indoor use, featuring four wheels for easy movement and a sturdy frame.

2. Safety Information

WARNING: Failure to follow these safety instructions may result in injury or damage to the device.

- Do not exceed the maximum weight capacity of 400 pounds (181 kg).
- Ensure proper assembly and use as instructed to avoid potential hazards.
- The device is designed for stability and safety but should be operated with caution.
- This lift chair is suitable for conscious users who can independently transfer onto the chair using hand strength.
- It is NOT recommended for individuals with paralysis, dementia, or those unable to self-support.
- This product is for indoor use only. Do not use outdoors.
- Not recommended for use on deep-pile or long-haired carpets, as this may impede wheel movement.
- Always lock the wheels before a user attempts to sit on or stand from the chair.





Important safety notice: This lift chair is suitable for conscious users who can independently transfer onto the chair using hand strength. It is not recommended for individuals with paralysis, dementia, or those unable to self-support.





This electric chair lift is designed exclusively for indoor use. It should not be used outdoors.

3. Package Contents

Verify that all components are present and undamaged upon opening the package:

- Electric Chair Lift (main unit)
- Armrests (2)
- Charging Adapter
- Lithium-Ion Battery (pre-installed or separate)
- Controller with Charging Line
- Safety Belt
- Assembly Tools (e.g., Allen wrench, screws)
- Instruction Manual




The package includes the charging adapter, battery, controller, and charging line necessary for operating the electric lift chair.

4. Setup and Assembly

The Jaioyios Electric Chair Lift features an enhanced foldable design for swift and effortless setup. Minimal assembly is required.

1. **Unfold the Chair:** Carefully unfold the main unit of the chair lift until it locks into an upright position.
2. **Attach Armrests:** Secure the two armrests to the designated slots on the chair frame using the provided screws and Allen wrench. Ensure they are firmly attached.
3. **Connect Controller and Battery:** Connect the controller's charging line to the battery. The battery may be pre-installed or require insertion into its compartment.
4. **Initial Charge:** Before first use, ensure the battery is fully charged. Refer to the 'Battery and Charging' section for details.



The chair lift folds down for convenient storage and transport, making it easy to move when not in use.



The armrests are designed to be foldable, facilitating easier transfers for the user and more compact storage of the chair.

Setup and Installation Video

This video demonstrates the assembly steps for the Jaioyios Electric Chair Lift, including attaching the armrests, and shows how to use its movable function to assist a person from the floor to a seated position.

5. Operating Instructions

The electric chair lift is designed for smooth and stable operation. Always ensure the user is properly secured before initiating any movement.

1. **Positioning the Chair:** Unfold the chair and move it to the desired location. The two front swivel wheels allow for 360-degree rotation, making it easy to maneuver.
2. **Locking the Wheels:** Before the user attempts to sit on the chair, press the brake levers on the front wheels to lock them, ensuring the chair remains stationary.
3. **Assisting the User:** Help the user transfer onto the chair. Ensure they are seated comfortably and securely.
4. **Secure with Safety Belt:** Fasten the buckle seat belt around the user to ensure stability during the lifting process.
5. **Operating the Lift:** Use the handheld controller to slowly and smoothly raise the chair. The controller allows for precise control over the lifting height.
6. **Unlocking Wheels for Movement:** Once the user is at a safe height and ready to be moved, unlock the front wheels by releasing the brake levers.
7. **Transferring the User:** Gently push the chair to transport the user to a bed, sofa, or other desired location.
8. **Lowering the Chair:** Use the controller to slowly lower the chair to the desired height for transfer or storage.



Step 1 of using the lift chair: Unfold the chair and position it near the person who needs assistance. Note: Not recommended for deep-pile or long-haired carpets.




Step 2: Lock the wheels for stability before the user sits on the chair. Step 3: Unlock the wheels to move the chair once the user is seated.



This image shows the electric lift chair in operation, smoothly raising a user from a seated position to a standing height, with a remote control for ease of use.



A buckle seat belt is included to ensure the user's comfort and stability during the lifting process.



The chair lift is equipped with two front swivel wheels that allow for smooth 360-degree rotation, enhancing maneuverability. These wheels can also be locked for stability.

6. Battery and Charging

The chair lift is powered by a high-capacity Lithium-Ion battery.

- **Battery Capacity:** The 2200mAh battery provides over 60 lifting cycles on a single full charge.
- **Charging Time:** A full charge typically takes 2-3 hours.
- **Charging Indicator:** The battery features a screen displaying the current charge level.

- **Charging Procedure:** Connect the charging adapter to the battery and plug it into a standard electrical outlet.



The chair lift is powered by a high-capacity 2200mAh Lithium-Ion battery, providing over 60 lifting cycles on a single charge. Charging time is 2-3 hours.

7. Maintenance

Regular maintenance ensures the longevity and safe operation of your chair lift.

- **Cleaning:** The PU seat and painted surfaces can be easily cleaned by wiping with a damp cloth. Avoid abrasive cleaners.

- **Battery Care:** Store the battery in a cool, dry place. If not in use for an extended period, charge the battery periodically to maintain its health.
- **Inspection:** Periodically inspect all moving parts, screws, and connections to ensure they are secure and free from damage.
- **Lubrication:** Consult a qualified technician for any necessary lubrication of mechanical components.

8. Troubleshooting

If you encounter issues with your Jaioyios Electric Chair Lift, refer to the following common solutions:

- **Chair Not Lifting:**
  - Check if the battery is fully charged. Recharge if necessary.
  - Ensure all cable connections (battery to controller, controller to chair) are secure.
  - Verify that the emergency stop button (if present) is not engaged.
- **Unstable Movement:**
  - Ensure all wheels are properly aligned and free of obstructions.
  - Check that the floor surface is even and suitable for wheeled movement (avoid deep carpets).
- **Armrests Loose:**
  - Tighten the screws securing the armrests using the provided Allen wrench.
- **Unusual Noises:**
  - Discontinue use immediately and inspect for any visible damage or loose parts. Contact customer support if the issue persists.

For issues not listed here, please contact Jaioyios customer support.

9. Specifications

Feature	Detail
Product Dimensions (Unfolded)	19"D x 17.7"W x 4"H (approximate, based on folded dimensions)
Product Dimensions (Folded)	19 x 17.7 x 4 inches
Item Weight	44.2 Pounds
Max Weight Capacity	400 Lbs
Model Number	AF1378-1
Battery Type	1 Lithium Ion battery (included)
Battery Capacity	2200mAh
Charging Time	2-3 hours
Lifting Cycles per Charge	Over 60

Material	Top-grade metal frame, PU seat
Wheels	4 wheels (2 front swivel with lock)
Color	Brown



This image illustrates the key dimensions of the chair lift, both unfolded and folded, highlighting its compact design for storage and its 400 lbs weight capacity.

10. Warranty and Support

For warranty information, please refer to the documentation included with your purchase or contact Jaioyios customer service directly. Our support team is available to assist with any questions or concerns regarding your electric chair lift.
Contact Information: [Insert Manufacturer Contact Information Here, if available]

© 2025 Jaioyios. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question