Rejuvacare RejuvaKnee 1.0

RejuvaKnee Triple Method Massager 1.0 Instruction Manual

Brand: Rejuvacare | Model: RejuvaKnee 1.0

1. Introduction

The RejuvaKnee Triple Method Massager 1.0 is designed to provide soothing comfort and relaxation for your knees. This device integrates gentle heat, adjustable compression, and therapeutic vibration to support overall knee well-being. It is suitable for everyday use to promote comfort, flexibility, and relaxation.



Image 1: RejuvaKnee Triple Method Massager 1.0 device and retail packaging.

2. Safety Information

Caution:

- Keep away from direct contact with skin.
- Do not use on broken skin or areas with inflammation.
- Avoid applying heat to sensitive areas.
- Consult a healthcare professional before use if you have any medical conditions or are pregnant.
- Keep the device away from water and moisture.
- Charge only with the provided charger.
- Do not use while sleeping or lying down.

3. Product Features

The RejuvaKnee 1.0 combines three methods for comprehensive knee support:

- **Soothing Heat:** Provides gentle warmth to comfort the knee area.
- **Vibration Massage:** Offers various levels of therapeutic vibration to help relieve stiffness and promote relaxation.
- **Customized Compression:** Adjustable straps allow for personalized pressure and stable support around the knee.



Image 2: Illustration of the Triple Therapy Technology (Heat, Air Compression, Vibration Massage).



Image 3: Effortless Control at Your Fingertips with touchscreen controls.

4. Setup

1. **Charge the Device:** Use the included USB charging cable to fully charge your massager before first use.
2. **Secure it around your Knee:** Place the massager over your knee, ensuring the control panel is positioned comfortably. Adjust the straps for a snug yet comfortable fit. The strong Velcro ensures secure placement.



Image 4: Charging the RejuvaKnee device.



Image 5: Securing the massager around the knee.

5. Operating Instructions

The RejuvaKnee 1.0 features a user-friendly one-button control interface:

- **Power On/Off:** Press and hold the power button (usually indicated by a circle with a vertical line) to turn the device on or off.
- **Adjust Heat:** Tap the '+' button to increase the heat level (indicated by a flame icon) or the '-' button to decrease it.
- **Adjust Massage Intensity:** Tap the wave icon button to cycle through different massage intensities or modes.



Image 6: Control Panel Diagram.

Recommended Use: Use the massager for 15 minutes daily for optimal comfort and relaxation.

6. Charging

The RejuvaKnee 1.0 is equipped with a rechargeable battery. To charge, connect the provided USB cable to the device's charging port and a suitable power source. Ensure the device is fully charged before use for best performance.

7. Maintenance

- **Cleaning:** Wipe the device with a soft, damp cloth. Do not immerse in water.
- **Storage:** Store the massager in a cool, dry place away from direct sunlight.

8. Troubleshooting

- **Device not turning on:** Ensure the battery is fully charged.
- **No heat/vibration:** Check that the device is powered on and settings are correctly adjusted.
- **Discomfort during use:** Re-adjust straps for a looser fit or reduce heat/massage intensity. Refer to safety guidelines.

9. Specifications

Specification	Detail
Product Dimensions	5.43 x 3.19 x 9.45 inches
Item Model Number	SF20250910B0FH2QZVQ8_B0FH2QZVQ8_1_HcYyD1
Batteries	1 Nonstandard Battery (included)
Use for	Knee
Power Source	Battery Powered
Material	Cotton
Item Weight	0.44 Grams
Brand	Rejuvacare

10. Official Product Videos

Watch these official videos for a visual guide on using your RejuvaKnee Triple Method Massager 1.0:

The Rejuvaknee Triple Method Massager

Your browser does not support the video tag.

Video 1: An official overview of the RejuvaKnee Triple Method Massager, demonstrating its features and benefits.

Knee Massager Smart with Heat

Your browser does not support the video tag.

Video 2: A demonstration of the smart knee massager with heat functionality.

11. Warranty & Support

For any questions, troubleshooting, or support regarding your RejuvaKnee Triple Method Massager 1.0, please contact Rejuvacare customer service. Refer to your product packaging or the official Rejuvacare website for contact details.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question