Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

KinseiBeauty T16

KinseiBeauty T16 Home Laser Hair Removal Device User Manual

1. Introduction

The KinseiBeauty T16 device is designed for gentle and effective hair removal in the comfort of your home. Utilizing advanced technology, it aims to provide smooth, hair-free skin and offers specialized functions for comprehensive skin care. This manual provides essential information for safe and optimal use of your device.

2. Safety Information

Please read all safety instructions before using the device. Failure to follow these instructions may result in injury or device damage.

- Always wear the provided eye protection during operation to prevent eye damage from intense light pulses.
- Do not use the device on tattooed areas, moles, birthmarks, or areas with skin conditions such as eczema, psoriasis, or open wounds.
- Consult a dermatologist before use if you have sensitive skin, a history of skin cancer, or are taking photosensitizing medications.
- Keep the device out of reach of children and pets.
- Do not use the device in wet environments or near water.
- Ensure the skin is clean, dry, and free of lotions, oils, or deodorants before treatment.
- Perform a patch test on a small skin area 24 hours before full treatment to check for adverse reactions.

3. Product Overview

The KinseiBeauty T16 device is designed for user comfort and efficiency. It features a lightweight, ergonomic body and comes with specialized lenses for different treatment purposes.



The KinseiBeauty T16 device features a sleek blue and gold design, engineered for comfortable handling during home use.

Components:

- Main Device Unit
- HR (Hair Removal) Lens
- SR (Skin Rejuvenation) Lens
- Power Adapter
- Protective Eyewear
- User Manual

4. Setup


Before your first treatment, follow these steps to set up your KinseiBeauty T16 device:

1. **Unpack:** Carefully remove all components from the packaging.
2. **Clean Skin:** Ensure the treatment area is clean, dry, and shaved. Avoid waxing or plucking, as the device targets the hair root.
3. **Connect Power:** Plug the power adapter into the device and then into an electrical outlet.
4. **Attach Lens:** Select the appropriate lens (HR for hair removal or SR for skin rejuvenation) and securely attach it to the device head.
5. **Wear Eye Protection:** Always put on the protective eyewear before powering on the device.

5. Operating Instructions

Follow these steps for effective and safe operation of your KinseiBeauty T16 device:

1. **Power On:** Press and hold the power button for a few seconds until the indicator lights illuminate.
2. **Select Intensity Level:** Short press the power button to cycle through the intensity levels. Start with the lowest level and gradually increase as your skin tolerates it.
3. **Apply to Skin:** Place the device window firmly against the skin in the desired treatment area. Ensure full contact.
4. **Flash:** Once the device detects full skin contact, the flash button will activate. Press the flash button to emit a light pulse.
5. **Move to Next Area:** Lift the device and move it to the next adjacent skin area. Avoid overlapping flashes excessively.
6. **Treatment Schedule:** For hair removal, typically use 1-2 times per week for the first 4-6 weeks, then reduce frequency as hair growth diminishes. For skin rejuvenation, follow the specific instructions for the SR lens.
7. **Power Off:** After completing your session, press and hold the power button to turn off the device.



This diagram illustrates how the laser targets hair follicles beneath the skin surface, inhibiting hair growth from the root. Removing hair from the roots can truly prevent hair growth, leaves no black spots, and prevents hair growth for a long time.

6. Treatment Modes (HR & SR)

The KinseiBeauty T16 comes with two specialized lenses to cater to different needs:

- **HR (Hair Removal) Lens:** This lens is designed to deliver light pulses that target melanin in the hair follicle, effectively inhibiting hair growth. It is suitable for various body areas including legs, arms, underarms, bikini line, and

face.

- **SR (Skin Rejuvenation) Lens:** This lens emits a different spectrum of light, optimized for skin rejuvenation. It can help improve skin texture, reduce the appearance of fine lines, and promote a more radiant complexion.

Always ensure you attach the correct lens for your desired treatment.

7. Maintenance

Proper maintenance ensures the longevity and optimal performance of your device:

- **Cleaning:** After each use, ensure the device is powered off and unplugged. Use a soft, dry cloth to wipe the treatment window and the device body. Do not use abrasive cleaners or immerse the device in water.
- **Lens Care:** Handle the lenses with care. If necessary, gently wipe the lens surface with a cotton swab moistened with alcohol.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and moisture. Keep it in its original packaging or a protective case when not in use.

8. Troubleshooting

If you encounter any issues with your KinseiBeauty T16 device, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Device does not power on	Not plugged in correctly; Power adapter fault	Check power connection; Try a different outlet; Contact support if adapter is faulty.
No flash emitted	Device not in full contact with skin; Low intensity level; Device malfunction	Ensure full skin contact; Increase intensity level; Contact support if issue persists.
Skin irritation or discomfort	Intensity level too high; Skin not properly prepared; Device used on sensitive areas	Reduce intensity; Ensure skin is clean and shaved; Avoid sensitive areas or use lowest setting. Apply a cool compress if irritation occurs.

If the problem persists after trying these solutions, please contact customer support.

9. Specifications

- **Model:** T16
- **Technology:** Intense Pulsed Light (IPL)
- **Modes:** HR (Hair Removal), SR (Skin Rejuvenation)
- **Lenses:** Two specialized lenses included
- **Cooling System:** Integrated Advanced Cooling
- **Flash Speed:** Ultra-Fast
- **Power Input:** AC 100-240V, 50/60Hz (Typical for home devices)
- **Color:** Blue (as per product description)

10. Warranty and Support

Your KinseiBeauty T16 device is backed by comprehensive support to ensure your satisfaction:

- **Warranty:** The device comes with a 2-year continuous technical support warranty.
- **Instruction Manual:** Includes a dermatologist's instruction manual for detailed guidance.
- **Follow-up Service:** Access to follow-up service for any questions or concerns regarding your device's performance or usage.

For warranty claims or technical assistance, please refer to the contact information provided in your product packaging or visit the official KinseiBeauty website.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0FMFQV6R4

More documents and user manuals: **https://manuals.plus/**