Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

CB CURABOX M156A

CB CURABOX Smart Automatic Pill Dispenser (Model M156A) Instruction Manual

Your guide to safe and effective medication management.

1. Introduction

The CB CURABOX Smart Automatic Pill Dispenser is designed to assist with organizing and dispensing medications, vitamins, and supplements. It features Wi-Fi connectivity for monitoring via a mobile application, audio/visual alerts, and a secure locking mechanism. This manual provides detailed instructions for setup, operation, and maintenance to ensure proper use.

2. Safety Information

- Keep this device out of reach of children.
- Do not use dosage rings as toys.
- Always use the lockable cover and manual key to prevent unauthorized access and accidental spills.
- Supervise children when they are near the pill dispenser.
- Ensure the device is placed on a stable, flat surface.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

3. Package Contents

Verify that all items are present in the package:

- CB CURABOX Smart Automatic Pill Dispenser (Model M156A)
- Locking Key
- USB Charging Cable
- Power Adapter
- Dosage Rings (9 included)
- Quick Start Guide
- Lithium Polymer Battery (pre-installed or separate)



Image: All components included with the CB CURABOX Smart Automatic Pill Dispenser.

4. Setup Guide

4.1. Unpacking and Initial Assembly

Carefully remove all components from the packaging. The dispenser has a transparent lid and an internal rotating tray.

4.2. Battery Installation and Power Connection

1. Locate the battery compartment on the underside of the dispenser.
2. Use the provided key to unlock the cover. Turn the key counter-clockwise to unlock.
3. Open the transparent lid and remove the internal medication tray.
4. Insert the Lithium Polymer battery into its compartment, ensuring correct polarity.
5. Connect the USB charging cable to the dispenser and the power adapter. Plug the adapter into a standard wall socket. The battery serves as a backup power source; keep the device charging during daily use.



Image: The dispenser's LCD screen displaying current date, time, Wi-Fi signal, battery status, and upcoming alarm information.

4.3. App Download and Device Pairing

1. Download the Tuya Smart Life app from your smartphone's app store.
2. Open the app and create an account or log in.
3. Ensure your phone is connected to a 2.4G Wi-Fi network (5G Wi-Fi is not supported for initial pairing).
4. In the app, use the auto-scan feature or manually add the device. Follow the on-screen prompts to connect the dispenser to your Wi-Fi network.



Image: The mobile app interface showing device connection and medication schedule monitoring.

4.4. Dosage Ring Placement

Select the appropriate dosage ring based on your medication schedule (e.g., 1x, 2x, 3x per day). Place the chosen ring onto the central mechanism of the dispenser.



Image: Customizable dosage rings for setting medication frequency.

4.5. Filling Medication Compartments

Fill each of the 28 high-capacity slots with the prescribed medication. Ensure to fill them in an anticlockwise direction, following your dosage schedule. The dispenser can accommodate multiple pills per slot.

4.6. Securing the Dispenser

Once filled, replace the internal tray and close the transparent lid. Use the key to lock the cover by turning it clockwise. This prevents accidental spills and unauthorized access.

5. Operating Instructions

5.1. Setting Medication Alarms

Use the mobile app to set up your medication schedule. Preset dosage templates provide multiple reminder options. You can adjust dispenser settings, set notification reminders, and alarms through the app.




Image: Setting timed reminders via the mobile app.

5.2. Receiving Medication Reminders

When it's time for medication, the dispenser will provide audio and visual alerts (light and voice). The mobile app will also send notifications to your phone.



Image: The dispenser's blue light indicator for medication reminders.

5.3. Dispensing Medication

To dispense medication, press the 'PRESS' button on the dispenser. The sliding medication access window will open, allowing you to retrieve your dose. The device will confirm the dose has been taken.




Image: Dispensing medication by pressing the 'PRESS' button.

5.4. Taking Medication in Advance

If you need to take a dose before its scheduled time, you can use the 'Take in advance' option within the mobile app. Confirm your action in the app, and the dispenser will release the next scheduled dose.

5.5. Reviewing Medication Records

The mobile app allows you to track your medication history. You can view monthly and daily records, indicating on-time doses (green), overdue doses (orange), and missed doses (red).



Image: Tracking medication records and compliance through the app.

5.6. Refilling Medication

When the dispenser indicates a refill is needed, unlock the cover, remove the tray, and replenish the medication compartments as described in the setup section.

5.7. Resetting the Dispenser Tray

After refilling the medication, it is crucial to reset the tray in the app. Navigate to the settings in the app and select the option to reset the tray. This synchronizes the dispenser with the new medication count.

5.8. Sharing Device Access

The app allows you to share device access with caregivers or family members. This enables them to monitor medication schedules and receive reminders, providing an additional layer of support.

Your browser does not support the video tag.

Video: A comprehensive overview of the CB CURABOX Smart Automatic Pill Dispenser's features, setup, and operation, including app integration and medication management.

6. Maintenance

6.1. Cleaning

To clean the dispenser, ensure it is unplugged and the battery is removed. Wipe the exterior with a soft, damp cloth. Do not submerge the device in water or use harsh chemicals. The internal tray can be gently wiped clean.

6.2. Battery Care

The built-in lithium battery provides backup power. For optimal performance, keep the dispenser connected to power during daily use. If storing the device for an extended period, ensure the battery is partially charged (around 50%) and store in a cool, dry place.

7. Troubleshooting

- **Device not powering on:** Ensure the battery is correctly installed and the device is connected to a working power outlet.
- **Wi-Fi connection issues:** Verify your router is broadcasting a 2.4G Wi-Fi signal. Ensure the dispenser is within range of the Wi-Fi router. Restart both the dispenser and your router if issues persist.
- **App not receiving notifications:** Check your phone's notification settings for the Smart Life app. Ensure the app has permission to send notifications.
- **Medication not dispensing:** Check if the dispenser is locked. Ensure the dosage ring is correctly placed and the compartments are not overfilled.
- **Incorrect time/date:** The device synchronizes with the app. Ensure the app is connected and the time settings are correct on your phone.

8. Specifications

- **Model Number:** M156A
- **Dimensions:** 10.04 x 9.84 x 2.95 inches
- **Weight:** 1.94 Pounds
- **Batteries:** 1 Lithium Polymer battery (included)
- **Material:** Plastic
- **Number of Compartments:** 28
- **Features:** Automatic, Built-In Alarm, Easy Open, Lockable, Rotatable, Wi-Fi Enabled

9. Support Information

For further assistance, please refer to the official CB CURABOX website or contact customer support through the Smart Life application. You may also find helpful resources and FAQs online.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0FN7B2Y6R

More documents and user manuals: **https://manuals.plus/**