Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

VEZURZO SR-07B

VEZURZO SR-07B 4-in-1 Electric Cellulite Massager & Body Sculpting Device User Manual

Model: SR-07B

1. Introduction

Thank you for choosing the VEZURZO SR-07B 4-in-1 Electric Cellulite Massager & Body Sculpting Device. This innovative device combines four advanced technologies—EMS (Electrical Muscle Stimulation), Red Light Therapy, Vacuum Suction, and Rolling Massage—to assist with body contouring, cellulite reduction, skin tightening, and improved circulation. Designed for convenience and effectiveness, the SR-07B offers a professional-grade experience in the comfort of your home.
Please read this manual thoroughly before first use to ensure safe and optimal operation. Keep this manual for future reference.

2. Safety Precautions

- **Application Gel:** Always apply a slimming gel, firming essential oil, or body lotion to the treatment area before use. This ensures smooth gliding of the device and maximizes comfort, preventing friction-related discomfort or skin irritation.
- **Intensity Levels:** Begin with the lowest intensity setting for both vacuum suction and EMS modes. Gradually increase the intensity as your skin adapts and you become comfortable with the sensation.
- **Treatment Duration:** Use the device for a maximum of 10 minutes per targeted area. The total session duration should not exceed 30 minutes. The device features an intelligent 10-minute auto-shutoff timer to prevent overuse.
- **Facial Treatment:** When using the small cup for facial treatment, exercise extreme caution and avoid the delicate eye area.
- **Cleaning:** Always clean the device thoroughly after each use to maintain hygiene, longevity, and effectiveness.
- **Medical Conditions:** Consult a healthcare professional before using this device if you are pregnant, have any underlying medical conditions, or are undergoing medical treatment.
- **Storage:** Store the device in a cool, dry place, out of reach of children.
- **Damage:** Do not use the device if it appears damaged or is not functioning correctly. Contact customer support for assistance.

3. Product Components



Figure 3.1: VEZURZO SR-07B device with both large and small cups.



Figure 3.2: The device comes with two interchangeable cups: a larger cup for body areas and a smaller cup for face and delicate areas.

The VEZURZO SR-07B device includes:

- Main Device Unit
- Large Suction Cup with Rolling Massagers
- Small Suction Cup with Rolling Massagers
- USB Charging Cable
- User Manual (this document)

4. Setup

4.1 Charging the Device

The device is powered by a 700mAh rechargeable battery. Before first use, fully charge the device. Connect the provided USB charging cable to the device and a standard USB power adapter (not included). The charging indicator light will show the charging status. A full

charge provides approximately 30 minutes of continuous use, sufficient for 3 treatment sessions.

4.2 Attaching the Suction Cup

1. Select the appropriate suction cup (large for body, small for face/arms) for your desired treatment area.
2. Align the chosen cup with the attachment point on the main device unit.
3. Gently push and twist the cup until it securely locks into place. Ensure it is firmly attached to prevent air leakage during operation.



Figure 4.1: The device with a cup attached, ready for use.

5. Operating Instructions

5.1 Preparing the Skin

Before turning on the device, apply a generous amount of slimming gel, firming essential oil, or body lotion to the skin area you intend to treat. This is crucial for the device to glide smoothly and for optimal contact and effectiveness of the EMS and vacuum functions.

5.2 Powering On/Off

- To power on the device, press and hold the power button (usually the central button) for approximately 2 seconds. The indicator lights will illuminate.
- To power off, press and hold the power button again for 2 seconds.



Figure 5.1: Device controls and indicator lights for power and mode selection.

5.3 Selecting Modes and Intensity

The SR-07B offers a 4-in-1 combination of EMS, Red Light Therapy, Vacuum Suction, and Rolling Massage. These functions can often be used simultaneously for enhanced results.

- **EMS (Electrical Muscle Stimulation) Modes:** The device features 5 EMS intensity modes. Press the EMS button (often labeled 'EMS' or indicated by a lightning bolt icon) to cycle through the intensity levels. Start with level 1 and increase as desired.
- **Vacuum Suction Levels:** There are 3 adjustable vacuum suction levels. Press the vacuum button (often labeled 'LIPO' or indicated by a suction cup icon) to select your preferred suction strength. Start with the lowest setting.
- **Red Light Therapy:** The red light therapy function activates automatically when the device is powered on and in use.
- **Rolling Massage:** The integrated rolling massagers work in conjunction with the vacuum suction to provide a deep tissue massage.




Figure 5.2: Overview of the 4-in-1 advanced technology features.

5.4 Treatment Application

Once the device is powered on and settings are selected, place the suction cup firmly against the prepared skin. Slowly glide the device over the treatment area. The vacuum suction will create a gentle pull on the skin, while the rollers and EMS work to stimulate the area.

- **Large Cup Areas:** Ideal for larger body areas such as the belly, thighs, buttocks, and love handles.
- **Small Cup Areas:** Suitable for smaller or more delicate areas like the face (avoiding eyes), double chin, arms, and neck.



Figure 5.3: Application of the device on the thigh area.



Figure 5.4: Application of the device on the face, using the smaller cup.

The device will automatically shut off after 10 minutes of continuous operation to ensure safe usage per area. For best results, consistent daily use for 2-3 minutes per area is recommended.

6. Maintenance

6.1 Cleaning the Device

Regular cleaning is essential for the hygiene and longevity of your VEZURZO SR-07B device.

1. Ensure the device is powered off and disconnected from the charging cable before cleaning.
2. Detach the suction cup from the main unit.
3. Wipe the exterior of the main unit with a soft, damp cloth. Do not submerge the main unit in water.
4. Clean the suction cups and rolling massagers with mild soap and water. Rinse thoroughly and allow them to air dry completely before reattaching or storing.
5. Avoid using abrasive cleaners, alcohol, or harsh chemicals, as these can damage the device.

6.2 Storage

Store the device in its original packaging or a clean, dry place away from direct sunlight and extreme temperatures. Ensure the device is fully charged before long-term storage to preserve battery health.

7. Troubleshooting

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or not charged.	Charge the device fully using the provided USB cable.
Weak or no vacuum suction.	Suction cup not properly attached; insufficient gel/oil on skin; low battery.	Ensure the cup is securely attached. Apply more gel/oil to the skin. Charge the device.
No EMS sensation.	EMS intensity too low; poor skin contact; insufficient conductive gel.	Increase EMS intensity. Ensure full contact with skin. Apply more conductive gel/lotion.

Problem	Possible Cause	Solution
Skin irritation or discomfort.	Insufficient gel/oil; intensity too high; prolonged use in one area.	Apply more gel/oil. Reduce intensity. Move the device continuously and limit use to 10 minutes per area.

If you encounter issues not listed here or if problems persist, please contact VEZURZO customer support for further assistance.

8. Specifications

Brand: VEZURZO

Model Name: SR-07B

Product Dimensions: 6 x 2.5 x 2.4 inches

Item Weight: 11 Ounces

Color: White

Material Type: Plastic

Power Source: Battery Powered (Rechargeable)

Battery: 1 A battery required (included), 700mAh

Cordless: Yes

Massage Therapy Technique: Cupping, Rolling

Special Features: Lightweight, Portable, Rechargeable, Auto-Shutoff

UPC: 762736230709

9. Warranty and Support

For warranty information, product support, or any inquiries regarding your VEZURZO SR-07B device, please refer to the contact information provided with your purchase or visit the official VEZURZO store on Amazon:

Visit the VEZURZO Store

Our customer service team is dedicated to assisting you with any questions or concerns.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B0FR8X6DSL

More documents and user manuals: **https://manuals.plus/**