Rejuvacare RejuvaFlow 2.0

Rejuvacare RejuvaFlow 3-in-1 Leg Massager 2.0 Instruction Manual

Model: RejuvaFlow 2.0

1. Important Safety Information

Please read all instructions carefully before using the RejuvaFlow 3-in-1 Leg Massager. Retain this manual for future reference.

- Consult your physician before use if you have any pre-existing medical conditions, including but not limited to, heart conditions, diabetes, pregnancy, or if you use a pacemaker.
- Do not use on swollen, inflamed, or injured areas of the skin.
- Keep the device away from water and high humidity.
- Do not use while sleeping.
- Keep out of reach of children.
- Use only the provided charging cable.
- Discontinue use immediately if you experience any discomfort or pain.

2. Package Contents

Ensure all items are present in your package:

- 2 x RejuvaFlow 3-in-1 Leg Massager Units
- 1 x USB Charging Cable
- 1 x User Manual



Figure 2.1: RejuvaFlow 3-in-1 Leg Massager 2.0 package contents.

3. Product Overview

The RejuvaFlow 3-in-1 Leg Massager 2.0 is designed to provide comfort and relaxation through a combination of three methods:

- **Massage:** Gentle stimulation to relax tired calves.
- **Compression:** Adjustable wrap for soothing support.
- **Heat:** Cozy warmth to help you unwind.



Figure 3.1: The RejuvaFlow 3-in-1 technology for everyday comfort.

The massager is cordless and portable, allowing for use in various settings such as home, office, or while traveling.



Figure 3.2: Enjoy the RejuvaFlow massager anytime, anywhere.



Figure 3.3: Versatile use of the RejuvaFlow massager in different environments.

4. Setup

1. **Charge the Units:** Before first use, fully charge both massager units using the provided USB cable. The display will indicate charging status.
2. **Wrap Around Calves:** Securely wrap each massager unit around your calves. Ensure a snug but comfortable fit using the adjustable straps. The control panel should face outwards for easy access.



Figure 4.1: Securely wrapping the massager around the calf.

5. Operating Instructions

1. **Power On:** Press and hold the power button on each unit for a few seconds until the display illuminates.
2. **Choose Mode:** Use the mode selection button to cycle through available massage, compression, and heat modes. Refer to the display for current settings.
3. **Adjust Intensity/Levels:** Use the '+' and '-' buttons to adjust the intensity of the massage or the level of heat/compression to your preference.
4. **Relax:** Enjoy 15-20 minutes of comfort. The device may have an automatic shut-off feature for safety.
5. **Power Off:** Press and hold the power button again to turn off the units.



Figure 5.1: Step-by-step guide to operating the RejuvaFlow massager.

The massager works by stimulating nerve endings, boosting circulation, and relieving tension in the legs.



Figure 5.2: Illustration of the massager's effect on leg anatomy.

6. Maintenance

- **Cleaning:** Wipe the massager units with a soft, damp cloth. Do not use abrasive cleaners or immerse in water.
- **Storage:** Store the massager in a cool, dry place, away from direct sunlight and extreme temperatures.
- **Battery Care:** To prolong battery life, charge the units regularly, even if not in frequent use.

7. Troubleshooting

- **Device not turning on:** Ensure the units are fully charged. Connect to the USB charger and try again.
- **No heat/vibration:** Check that the desired mode is selected and the intensity level is adjusted appropriately.
- **Discomfort during use:** Reduce the intensity level or adjust the tightness of the wrap. If discomfort persists, discontinue use.

8. Specifications

Feature	Detail
Product Dimensions	5.43"L x 3.19"W x 9.45"H
Brand	Rejuvacare
Color	Grey
Material	Cotton
Special Feature	Cordless
Manufacturer	Rejuvacare
ASIN	B0FX5N99G8

9. Warranty and Support

For warranty information or customer support, please refer to the contact details provided on the product packaging or visit the official Rejuvacare website. Keep your purchase receipt for warranty claims.

© 2025 Rejuvacare. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for

conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question