Rejuvacare 2026

Rejuvacare RejuvaKnee Triple Method Massager Instruction Manual

Model: 2026

Introduction

Thank you for choosing the Rejuvacare RejuvaKnee Triple Method Massager. This device is designed to provide comfort and relaxation for your knees through a combination of heat and massage therapy. Please read this manual thoroughly before use to ensure safe and effective operation.



Image: The RejuvaKnee Triple Method Massager, a grey knee wrap with a black control unit on top.

Important Safety Information

Always follow basic safety precautions when using electrical appliances.

- Do not use if you have any medical conditions that contraindicate massage or heat therapy. Consult a physician before use if you have pre-existing conditions, are pregnant, or have a pacemaker.
- Keep out of reach of children.
- Do not immerse the device in water or use it in wet environments.
- Avoid using the device on swollen, inflamed, or injured areas of the skin.
- Do not use for more than 30 minutes per session.
- Ensure the charging cable is not damaged before use.
- Use only the provided charging accessories.

Package Contents

- 1 x RejuvaKnee Triple Method Massager (Knee Wrap with integrated control unit)
- 1 x Charging Cable
- 1 x User Manual (This document)

Product Overview

The RejuvaKnee Triple Method Massager is a wearable knee wrap designed for targeted comfort. It features an intuitive control panel for adjusting settings.



Image: A close-up view of the RejuvaKnee control panel, showing the power button, wave icon for massage intensity, and plus/minus buttons for heat adjustment.

Control Panel Functions:

- **Power Button (Circle with vertical line):** Press and hold to turn the device on or off.
- **Massage Intensity Button (Wavy lines):** Tap to cycle through different massage intensity levels.
- **Heat Increase Button (+):** Tap to increase the heat level.
- **Heat Decrease Button (-):** Tap to decrease the heat level.

Setup

1. **Charge the Device:** Before first use, fully charge the RejuvaKnee massager. Connect the charging cable to the device and a suitable power source. The indicator light will show charging status.
2. **Adjust the Wrap:** Place the RejuvaKnee wrap around your knee. Use the adjustable straps to secure it comfortably but firmly. Ensure the control unit is positioned on the front of your knee for easy access.



Image: A person adjusting the RejuvaKnee massager on their knee, demonstrating the adjustable wrap-around fit.

Operating Instructions

1. **Power On:** Press and hold the Power button on the control unit for a few seconds until the device turns on.
2. **Select Massage Intensity:** Tap the Massage Intensity button (wavy lines) to choose your desired vibration level.
3. **Adjust Heat Level:** Use the '+' and '-' buttons to increase or decrease the heat intensity.
4. **Automatic Timer:** The device is equipped with a 30-minute auto-timer. It will automatically shut off after 30 minutes of continuous use. This is to prevent overuse and ensure safety.
5. **Power Off:** To turn off the device manually, press and hold the Power button again until it shuts down.

Your browser does not support the video tag.

Video: Official product video demonstrating the features and use of the RejuvaKnee Triple Method Massager.



Image: A person using the RejuvaKnee massager on their knee while relaxing, highlighting its ease of use during daily activities.

Charging

The RejuvaKnee massager is cordless and rechargeable.

- Connect the provided USB charging cable to the charging port on the control unit.
- Plug the other end of the cable into a standard USB power adapter (not included) or a computer USB port.
- The indicator light on the control unit will show the charging status. Refer to the device's display for battery level.
- A full charge typically takes approximately 2-3 hours.
- Do not operate the device while it is charging.



Image: The RejuvaKnee massager showing its wireless design and charging port, emphasizing its portability.

Maintenance and Care

- **Cleaning:** Wipe the surface of the massager with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the device in a cool, dry place away from direct sunlight and extreme temperatures.
- **Battery Care:** To prolong battery life, avoid fully discharging the battery frequently. Charge it regularly, even if not in use for extended periods.

Troubleshooting

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or device is off.	Charge the device fully. Press and hold the power button to turn on.

Problem	Possible Cause	Solution
No heat or massage.	Settings not activated or device malfunction.	Ensure heat and massage functions are selected using their respective buttons. If problem persists, contact customer support.
Device stops unexpectedly.	Auto-timer activated or low battery.	The device has a 30-minute auto-timer. Recharge the device if the battery is low.
Charging indicator not working.	Cable not connected properly or faulty cable/charger.	Ensure cable is securely connected. Try a different USB port or charger.

Specifications

- **Model Name:** RejuvaKnee
- **Model Number:** 2026
- **Product Dimensions:** 5.43"L x 3.19"W x 5.45"H
- **Item Weight:** 0.44 Grams (0.02 ounces)
- **Power Source:** Battery Powered (1 Nonstandard Battery included)
- **Target Use Body Part:** Knee
- **Product Benefits:** Pain Relief, Muscle Therapy
- **Special Features:** Triple Method Massager (Kneading, Percussion, Rolling), Cordless
- **Material Type:** Cotton
- **Color:** Grey
- **Manufacturer:** Rejuvacare

Warranty and Support

Rejuvacare offers a 90-day promise on its products. For warranty claims, technical support, or any questions regarding your RejuvaKnee Triple Method Massager, please contact Rejuvacare customer service.
Please visit the official Rejuvacare website or refer to your purchase documentation for contact details.

Disposal Information

This product contains a battery and electronic components. Do not dispose of it with general household waste. Please follow local regulations for the proper disposal of electronic waste and batteries. Contact your local waste management authority for information on recycling options.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question