Rejuvacare Foot Renew

Rejuvacare Foot Renew Triple Method Massager Instruction Manual

Model: Foot Renew | Brand: Rejuvacare

Introduction

The Rejuvacare Foot Renew Triple Method Massager is designed to provide a soothing experience for your feet and ankles. This device combines heat, vibration, and compression to offer personalized comfort and support for daily relaxation at home. Please read this manual thoroughly before using the device to ensure proper operation and safety.

Product Features

- **Foot & Ankle Comfort:** Designed to provide a soothing foot and ankle massage experience with adjustable settings.
- **Adjustable Modes:** Select your preferred mode and level to match your routine, created for unwinding after long days.
- **Comfortable & Easy Setup:** Simple to position and use, fitting easily into your daily routine.
- **Suitable for Daily Relaxation:** Ideal for relaxing while watching TV, working from home, or winding down before bed.
- **Compact Design:** Includes 1 Foot Renew foot & ankle massager, easy to store.



Image: The Rejuvacare Foot Renew massager, showcasing its compact design and overall appearance.



Image: Diagram illustrating the three methods of therapy: mild heat, compression, and vibration massage, targeting foot discomfort.

Setup

Follow these steps to prepare your Rejuvacare Foot Renew massager for first use:

1. **Charge the Device:** Connect the massager to a power source using the included USB-C cable. The display will indicate charging status.
2. **Position on Foot/Ankle:** Place the massager on your foot or ankle as needed, ensuring a comfortable and secure fit. The adjustable strap allows for various foot sizes.

Image: Visual guide demonstrating the simple steps to charge the device, put it on the foot, and operate the controls for heating and vibration.

Operating Instructions

To operate your Rejuvacare Foot Renew massager, use the control panel located on the device:

- **Power On/Off:** Press and hold the Power button ( ) for 2 seconds to turn the device on or off. The default temperature is 50°C.
- **Adjust Heating:** Press the Heating button ( ) to cycle through available temperature settings.
- **Adjust Vibration:** Press the Vibration button (≋) to adjust the vibration mode and frequency.

Your browser does not support the video tag.

Video: An official demonstration by Rejuvacare US showing how to wear and operate the Foot Renew Wireless Electric Foot & Calf Massager.

Care and Maintenance

To ensure the longevity and proper functioning of your Rejuvacare Foot Renew massager, follow these care guidelines:

- **Cleaning:** Wipe the device with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- **Storage:** Store the massager in a cool, dry place away from direct sunlight and extreme temperatures.
- **Battery Care:** Charge the device regularly, even if not in frequent use, to maintain battery health.

Troubleshooting

If you encounter any issues with your Rejuvacare Foot Renew massager, refer to the following common solutions:

- **Device Not Turning On:** Ensure the device is fully charged. Press and hold the power button for at least 2 seconds.
- **No Heat/Vibration:** Check if the respective function is activated on the control panel. Ensure the device has sufficient battery charge.
- **Uncomfortable Fit:** Adjust the strap to find a comfortable and secure fit around your foot or ankle.

If problems persist, please contact Rejuvacare customer support.

Specifications

Feature	Detail
Product Dimensions	3 x 9 x 5 inches
Weight	11.36 ounces
Manufacturer	Rejuvacare
Target Use Body Part	Foot
Special Features	Lightweight, Portable
Unit Count	1 Count
Color	Black
Material Type	Polyester

Safety Information

Please observe the following safety precautions when using the Rejuvacare Foot Renew massager:

- Do not use the device if you have any open wounds, skin irritations, or severe circulatory problems in the area of application.
- Consult a healthcare professional before use if you have pre-existing medical conditions, are pregnant, or have a pacemaker.
- Keep the device away from water and moisture.
- Do not use while sleeping.
- Keep out of reach of children.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This device is for relaxation and comfort support.

Warranty and Support

Rejuvacare is committed to providing quality products. For any questions regarding your Foot Renew massager, including warranty information or technical support, please refer to the contact information provided with your purchase or visit the official Rejuvacare website.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Answered questions about this manual

What settings should I use for heat and vibration on the Rejuvacare Foot Renew massager?

The Rejuvacare Foot Renew massager features adjustable settings for both heat and vibration to customize your experience. Here is how to use them:Power On: Press and hold the Power button ( ) for 2 seconds....

Read full answer